



PRIDE JAZZY SELECT ELITE HD

~~\$5,195.99~~ \$4,095.99

The Jazzy® Elite HD features two-motor, front-wheel drive technology and heavy-duty construction.

The Jazzy Elite HD features 14" drive wheels for enhanced outdoor performance combined with a front-wheel drive design to provide tight turns for optimal indoor maneuverability in small spaces.

Large front wheels add extra absorption for superb climbing capabilities and combine with front antitip wheels for handling various terrains with ease.

- Enhanced 204kg (450 lbs.) weight capacity
- Front-wheel drive design for excellent maneuverability
- 14" knobby tires and optimum weight distribution provide superior traction
 - Larger foot platform
 - 60 amp controller



SKU: S10463

Category: [Power Chairs](#)

ADDITIONAL INFORMATION

Brand	Pride
Max Weight capacity (Max SWL)	204 kg (450 lbs)
Maximum speed	Up to 6.4 km/h (4.0 mph)
Ground clearance	7.6 cm at bottom of base (3 inches)
Turning radius	61 cm with front riggings (24 inches)
High-Back Seat (max. size)	60cm W x 60cm D 24 inches W x 24 inches D
Standard seat size	50.8 x 45.7 – 50.8cm 20 x 18-20 inches
Overall length	90.2 cm without front riggings (35.5 inches)
Overall width	70 cm (27.5 inches)
Drive wheels	35.6cm (14 inches foam filled)
Warranty	2 years