




**SKU:** S2446

## **PRIDE LC-101 ELECTRIC RECLINER LIFT CHAIR - LEATHER**



**\$2,299.99**

A Pride LC-101 is a three position, single motor lift chair with classic style and exceptional quality and performance. Pride's LC-101 delivers superb value, excellent quality and luxurious comfort. The Pride LC-101 lift chair features quiet and smooth operation and stylish design to add to any home décor.



Scan QR code to  
view product online

VARIATIONS

Image	SKU	Price	Description	Leather
	S2446	\$2,299.99		Beige
	S2447	\$2,299.99		Brown
	S6913	\$2,299.99		Black

## PRODUCT DESCRIPTION

The Pride LC-101 is a three-position, single-motor lift chair with a classic design and high quality and performance. A handy integrated emergency battery back-up system ensures mobility in the event of a power outage.

### **Single vs twin vs quad motor lift chairs - what's the difference?**

All ILS lift chairs have a lifting action, to support the user to stand from a seated position.

- Single motor chairs move the back rest and leg rest at the same time.
- Twin motor chairs allow the backrest and leg rest to move independently of each other, typically allowing for a greater degree of leg elevation and greater degree of recline compared to a standard single motor chair.
- Quad motor chairs allow independent movement of headrest, lumbar, backrest and footrest allowing for the greatest level of re-positioning of all chairs.

## ADDITIONAL INFORMATION

<b>Brand</b>	Pride
<b>Seat Height</b>	470 mm
<b>Seat Depth</b>	510 mm
<b>Seat Width</b>	510 mm
<b>Overall Width</b>	900 mm
<b>Overall Height</b>	1040 mm
<b>Max SWL</b>	147 kg