



**SKU:**

## **PRIDE MOBILITY LC-101 EURO LEATHER LIFT CHAIR (SINGLE MOTOR)**





**\$2,699.99**

The Pride Mobility LC-101 Euro Leather Lift Chair (Single Motor) provides practical support in a sophisticated finish. Its three-position mechanism and rounded European armrests offer both functionality and design appeal.



Scan QR code to  
view product online

VARIATIONS

Image	SKU	Price	Description	Leather
	S2418	\$2,699.99		Burgundy
	S2446	\$2,699.99		Beige
	S2447	\$2,699.99		Brown
	S6913	\$2,699.99		Black

## PRODUCT DESCRIPTION

The Pride Mobility LC-101 Euro Leather Lift Chair (Single Motor) is a stylish, single motor model that transitions between seated, reading and lift positions with ease.

The chair's European-style rounded armrests and high-quality euro leather finish create a refined appearance, while its mechanism supports safe transfers. There are four modern colour options to choose from.

Equipped with an intuitive hand control featuring USB charging and battery backup, it's a dependable and attractive lift solution.

### Who it's for

- For clients who want a leather lift chair with everyday recline and lift functions.
- For therapists needing a functional yet elegant option for residential environments.

### Key features

- **Three-position single motor design** – supports seated, reading and standing positions
- **Elegant euro leather upholstery** – stylish, soft and easy to maintain
- **Rounded armrests** – classic European look
- **Durable internal frame** – built with quality wood and steel
- **Easy-to-use hand control with USB** – easy to operate and charge devices
- **Battery backup** – ensures consistent functionality
- **Side pockets** – for convenient personal storage

## ADDITIONAL INFORMATION

<b>Brand</b>	Pride Mobility
<b>Seat Height</b>	470 mm
<b>Seat Depth</b>	510 mm
<b>Seat Width</b>	510 mm
<b>Overall Width</b>	900 mm
<b>Overall Height</b>	1040 mm
<b>Max SWL</b>	147 kg