




**SKU:**

## **PRIDE LC-358XL BARIATRIC ELECTRIC RECLINER LIFT CHAIR**



**\$3,699.99**

Comfort and style complement one another perfectly in the LC-358XL from Pride® Lift Chairs. This lift chair has a superior design that provides exceptional strength, style, and durability that will last for years. The LC-358 XL 3-position chaise lift chair from Pride Mobility offers ample leg room, courtesy of a padded footrest extension. Make this one-of-a-kind power lift recliner part of your home! The LC-358 XL also benefits from a heavy-duty 227 kg weight capacity. Superb comfort, smooth and quiet performance all at a great value! This lift chair comes in two variations: Latte (Fabric) and Burgundy (Vinyl).



Scan QR code to  
view product online

VARIATIONS

Image	SKU	Price	Description	Colours
	S9829	\$3,699.99		Latte (Fabric)
	S9828	\$3,699.99		Burgundy (Vinyl)

## PRODUCT DESCRIPTION

Classic style meets exquisite comfort. The Single Motor LC-358 XL 3-position chaise lift chair from Pride Mobility offers ample leg room courtesy of a padded footrest extension. Make this one-of-a-kind power lift recliner part of your home! The LC-358 XL also benefits from a heavy-duty 227 kg weight capacity.

- The 3-Position Full Recline puts you in a position ideal for napping or relaxing
- Exceptional strength, durability, and superior design
- Heavy Duty (SWL 227 kg)

### Single vs twin vs quad motor lift chairs - what's the difference?

All ILS lift chairs have a lifting action, to support the user to stand from a seated position.

- Single motor chairs move the back rest and leg rest at the same time.
- Twin motor chairs allow the backrest and leg rest to move independently of each other, typically allowing for a greater degree of leg elevation and greater degree of recline compared to a standard single motor chair.
- Quad motor chairs allow independent movement of headrest, lumbar, backrest and footrest allowing for the greatest level of re-positioning of all chairs.

## ADDITIONAL INFORMATION

<b>Weight Capacity / Max SWL</b>	227 kg
<b>Back Width</b>	91 cm
<b>Seat To Floor</b>	51 cm
<b>Seat Depth</b>	56 cm
<b>Seat Width</b>	66 cm
<b>Top Of Back To Seat</b>	74 cm
<b>Integrated Battery Backup</b>	Yes
<b>Fabric Selection Options</b>	Latte (Fabric) & Burgandy (Vinyl)
<b>Motor Type</b>	Single Motor