



SKU:

PRIDE MOBILITY LC-358XL LIFT CHAIR (SINGLE MOTOR)

\$3,899.99

The Pride Mobility LC-358XL Lift Chair (Single Motor) is a robust seating solution supporting up to 227 kg. Featuring a heavy-duty frame, it combines supportive foam cushioning with simple single-motor operation to assist with repositioning and transfers.



Scan QR code to
view product online

VARIATIONS

Image	SKU	Price	Description	Colours
	S9829	\$3,899.99		Latte (Fabric)
	S9828	\$3,899.99		Burgundy (Vinyl)

PRODUCT DESCRIPTION

The Pride LC-358XL Lift Chair (Single Motor) is a bariatric lift chair designed for individuals who need a higher weight capacity and reliable seating support.

The single motor mechanism allows the backrest and legrest to move together, creating a smooth and quiet transition between reclined and upright positions. This supports safer transfers and enhances daily comfort.

Built on a heavy-duty hardwood frame, the chair features high-density foam cushioning to enhance comfort and provide pressure relief.

The chair also includes an easy-to-use hand control with USB charging.

Available in Latte Suede fabric or Burgundy Vinyl upholstery.

Who it's for

- For clients who need bariatric seating with smooth recline and lift support.
- For therapists wanting a durable lift chair for clients up to 227 kg.

Key features

- **Single motor operation** – reclines and lifts for smooth, simultaneous positioning
- **Heavy-duty hardwood frame** – supports individuals up to 227 kg
- **High-density foam cushioning** – enhances long-term comfort and pressure relief
- **Quiet, smooth mechanism** – for stress-free use
- **Easy-to-use hand control with USB** – keeps personal devices powered
- **Battery backup** – maintains function during power interruptions
- **Side pockets** – for convenient personal storage

ADDITIONAL INFORMATION

Brand	Pride Mobility
Weight Capacity / Max SWL	227 kg
Back Width	91 cm
Seat To Floor	51 cm
Seat Depth	56 cm
Seat Width	66 cm
Top Of Back To Seat	74 cm
Integrated Battery Backup	Yes
Fabric Selection Options	Latte (Fabric) & Burgandy (Vinyl)