




SKU:

PRIDE LC-358XL BARIATRIC ELECTRIC RECLINER LIFT CHAIR

\$3,899.99

Comfort and style complement one another perfectly in the LC-358XL from Pride® Lift Chairs. This lift chair has a superior design that provides exceptional strength, style, and durability that will last for years. The LC-358 XL 3-position chaise lift chair from Pride Mobility offers ample leg room, courtesy of a padded footrest extension. Make this one-of-a-kind power lift recliner part of your home! The LC-358 XL also benefits from a heavy-duty 227 kg weight capacity. Superb comfort, smooth and quiet performance all at a great value! This lift chair comes in two variations: Latte (Fabric) and Burgundy (Vinyl).



Scan QR code to
view product online

VARIATIONS

Image	SKU	Price	Description	Colours
	S9829	\$3,899.99		Latte (Fabric)
	S9828	\$3,899.99		Burgundy (Vinyl)

PRODUCT DESCRIPTION

Classic style meets exquisite comfort. The Single Motor LC-358 XL 3-position chaise lift chair from Pride Mobility offers ample leg room courtesy of a padded footrest extension. Make this one-of-a-kind power lift recliner part of your home! The LC-358 XL also benefits from a heavy-duty 227 kg weight capacity.

- The 3-Position Full Recline puts you in a position ideal for napping or relaxing
- Exceptional strength, durability, and superior design
- Heavy Duty (SWL 227 kg)

Single vs twin vs quad motor lift chairs - what's the difference?

All ILS lift chairs have a lifting action, to support the user to stand from a seated position.

- Single motor chairs move the back rest and leg rest at the same time.
- Twin motor chairs allow the backrest and leg rest to move independently of each other, typically allowing for a greater degree of leg elevation and greater degree of recline compared to a standard single motor chair.
- Quad motor chairs allow independent movement of headrest, lumbar, backrest and footrest allowing for the greatest level of re-positioning of all chairs.

ADDITIONAL INFORMATION

Weight Capacity / Max SWL	227 kg
Back Width	91 cm
Seat To Floor	51 cm
Seat Depth	56 cm
Seat Width	66 cm
Top Of Back To Seat	74 cm
Integrated Battery Backup	Yes
Fabric Selection Options	Latte (Fabric) & Burgandy (Vinyl)
Motor Type	Single Motor
Brand	Pride