



**SKU:**

## **PRIDE MOBILITY LC-358XL LIFT CHAIR (SINGLE MOTOR)**

**\$3,899.99**

The Pride Mobility LC-358XL Lift Chair (Single Motor) is a robust seating solution supporting up to 227 kg. Featuring a heavy-duty frame, it combines supportive foam cushioning with simple single-motor operation to assist with repositioning and transfers.



Scan QR code to  
view product online

## VARIATIONS

| Image  | SKU   | Price      | Description | Colours          |
|--|-------|------------|-------------|------------------|
|   | S9829 | \$3,899.99 |             | Latte (Fabric)   |
|  | S9828 | \$3,899.99 |             | Burgundy (Vinyl) |

## PRODUCT DESCRIPTION

The Pride LC-358XL Lift Chair (Single Motor) is a bariatric lift chair designed for individuals who need a higher weight capacity and reliable seating support.

The single motor mechanism allows the backrest and legrest to move together, creating a smooth and quiet transition between reclined and upright positions. This supports safer transfers and enhances daily comfort.

Built on a heavy-duty hardwood frame, the chair features high-density foam cushioning to enhance comfort and provide pressure relief.

The chair also includes an easy-to-use hand control with USB charging.

Available in Latte Suede fabric or Burgundy Vinyl upholstery.

### Who it's for

- For clients who need bariatric seating with smooth recline and lift support.
- For therapists wanting a durable lift chair for clients up to 227 kg.

### Key features

- **Single motor operation** – reclines and lifts for smooth, simultaneous positioning
- **Heavy-duty hardwood frame** – supports individuals up to 227 kg
- **High-density foam cushioning** – enhances long-term comfort and pressure relief
- **Quiet, smooth mechanism** – for stress-free use
- **Easy-to-use hand control with USB** – keeps personal devices powered
- **Battery backup** – maintains function during power interruptions
- **Side pockets** – for convenient personal storage

## ADDITIONAL INFORMATION

|                                  |                                   |
|----------------------------------|-----------------------------------|
| <b>Brand</b>                     | Pride Mobility                    |
| <b>Weight Capacity / Max SWL</b> | 227 kg                            |
| <b>Back Width</b>                | 91 cm                             |
| <b>Seat To Floor</b>             | 51 cm                             |
| <b>Seat Depth</b>                | 56 cm                             |
| <b>Seat Width</b>                | 66 cm                             |
| <b>Top Of Back To Seat</b>       | 74 cm                             |
| <b>Integrated Battery Backup</b> | Yes                               |
| <b>Fabric Selection Options</b>  | Latte (Fabric) & Burgandy (Vinyl) |