




SKU:

PRIDE LC 107 TWIN MOTOR LIFT CHAIR




\$1,999.99

The LC-107 is one of Pride's most versatile lift chairs. The chair is infinite position allowing for a variety of comfort, from watching TV to sleeping. The LC-107 is available at a low cost, making it an exceptional value.



Scan QR code to
view product online

VARIATIONS

Image	SKU	Price	Description	Fabric
	S10077	\$1,999.99		Burgundy
	S10079	\$1,999.99		Ocean
	S10078	\$1,999.99		Pecan

PRODUCT DESCRIPTION

Look no farther than Pride's infinite-position LC-107 power lift recliner for a versatile lift chair that will provide therapeutic support. This beautiful model is ideal for unwinding in front of the television and offers excellent value for money.

Single vs twin vs quad motor lift chairs - what's the difference?

All ILS lift chairs have a lifting action, to support the user to stand from a seated position.

- Single motor chairs move the back rest and leg rest at the same time.
- Twin motor chairs allow the backrest and leg rest to move independently of each other, typically allowing for a greater degree of leg elevation and greater degree of recline compared to a standard single motor chair.
- Quad motor chairs allow independent movement of headrest, lumbar, backrest and footrest allowing for the greatest level of re-positioning of all chairs.

ADDITIONAL INFORMATION

Brand	Pride
Seat To Floor	510 mm
Seat Depth	530 mm
Seat Width	510 mm
Top of Back to Seat	700 mm
Max SWL	147 kg
Motor	Twin