

SKU:

## PRIDE MOBILITY LC-107 LIFT CHAIR (TWIN MOTOR)

**\$2,299.99**

The Pride Mobility LC-107 Lift Chair (Twin Motor) offers independent control of the backrest and legrest for a customisable seating experience. With a medium seat size and classic style, it's a dependable choice for everyday lift support.



Scan QR code to  
view product online

### VARIATIONS

SKU	Price	Description	Colours
S10077	\$2,299.99	The Pride Mobility LC-107 Lift Chair (Twin Motor) offers independent control of the backrest and legrest for a customisable seating experience. With a medium seat size and classic style, it's a dependable choice for everyday lift support.	Burgundy
S10079	\$2,299.99	The Pride Mobility LC-107 Lift Chair (Twin Motor) offers independent control of the backrest and legrest for a customisable seating experience. With a medium seat size and classic style, it's a dependable choice for everyday lift support.	Ocean
S10078	\$2,299.99	The Pride Mobility LC-107 Lift Chair (Twin Motor) offers independent control of the backrest and legrest for a customisable seating experience. With a medium seat size and classic style, it's a dependable choice for everyday lift support.	Pecan

## PRODUCT DESCRIPTION

The Pride Mobility LC-107 Lift Chair (Twin Motor) is designed for people who want more control over their seating position.

Its dual-motor system allows independent adjustment of the backrest and legrest, helping individuals find their most comfortable position with ease.

This medium-sized chair features high-density foam, extended footrest, and an easy-to-use USB hand control. Available in three modern colour options.

### Who it's for

- For clients who need a flexible and supportive lift chair.
- For therapists wanting dual-motor positioning and simple functionality in a medium-sized option.

### Key features

- **Twin motor operation** - enables independent adjustment of backrest and legrest
  - **Medium-size seat** - fits a wide range of individuals
  - **Extended footrest** - supports reclined positions and taller individuals
  - **Durable frame and steel mechanism** - built for reliability
- **User-friendly hand control with USB** - easy to operate and charge devices
  - **Battery backup** - ensures function during power outages
  - **Side pockets** - for convenient storage

## ADDITIONAL INFORMATION

<b>Safe Working Load (SWL)</b>	147 kg
<b>Overall Width</b>	109 cm
<b>Overall Depth</b>	90 cm
<b>Overall Height</b>	108 cm
<b>Seat Width</b>	51 cm
<b>Seat Depth</b>	53 cm
<b>Seat to Floor Height</b>	51 cm
<b>Number of Motors</b>	Twin
<b>Brand</b>	Pride Mobility

[Pride Mobility Brochure](#)

[Pride Mobility LC-107 User Manual](#)

[Pride Mobility Range Comparison Chart](#)

### **Talk to us on the phone**

Our experienced phone sales team can assist you with any product enquiries and help you make a purchase directly on the phone. We are available anytime between 9am and 5pm, Monday to Friday.

[1300 008 267](tel:1300008267)

**[SINGLE VS TWIN VS QUAD MOTOR LIFT CHAIRS - WHAT'S THE DIFFERENCE?](#)**