




**SKU:**

## **PRIDE LC 107 TWIN MOTOR LIFT CHAIR**




**\$1,999.99**

The LC-107 is one of Pride's most versatile lift chairs. The chair is infinite position allowing for a variety of comfort, from watching TV to sleeping. The LC-107 is available at a low cost, making it an exceptional value.



Scan QR code to  
view product online

### VARIATIONS

Image	SKU	Price	Description	Fabric
	S10077	\$1,999.99		Burgundy
	S10079	\$1,999.99		Ocean
	S10078	\$1,999.99		Pecan

## PRODUCT DESCRIPTION

Look no farther than Pride's infinite-position LC-107 power lift recliner for a versatile lift chair that will provide therapeutic support. This beautiful model is ideal for unwinding in front of the television and offers excellent value for money.

### **Single vs twin vs quad motor lift chairs - what's the difference?**

All ILS lift chairs have a lifting action, to support the user to stand from a seated position.

- Single motor chairs move the back rest and leg rest at the same time.
- Twin motor chairs allow the backrest and leg rest to move independently of each other, typically allowing for a greater degree of leg elevation and greater degree of recline compared to a standard single motor chair.
- Quad motor chairs allow independent movement of headrest, lumbar, backrest and footrest allowing for the greatest level of re-positioning of all chairs.

## **ADDITIONAL INFORMATION**

<b>Brand</b>	Pride
<b>Seat To Floor</b>	510 mm
<b>Seat Depth</b>	530 mm
<b>Seat Width</b>	510 mm
<b>Top of Back to Seat</b>	700 mm
<b>Max SWL</b>	147 kg
<b>Motor</b>	Twin