




SKU: S2224

PRIDE PATHRIDER 10 DELUXE MOBILITY SCOOTER

From: \$3,649.99

The sleek, sporty Pathrider 10 Deluxe scooter delivers smooth performance and all-new features. With a full complement of advanced standard features including, rugged suspension, a high-visibility automatic front and rear lighting system, feather-touch disassembly, an LED kerb light, and a wraparound easy-drive tiller, the Pathrider 10 Deluxe delivers greater value than ever before.



Scan QR code to
view product online

PRODUCT DESCRIPTION

Key Features

- 159 kg weight capacity
- 157.5 cm turning radius
- Up to 8.5 km/h maximum speed
- Up to 25 km per charge
- Exclusive one-hand feather touch disassembly provides the ultimate quick transport system
- Easy-drive tiller with wraparound handles
- Back-lit battery gauge
- Front and rear suspension
- Lightweight pinch less hinge seat, eight position rotation, sliders, and viscoelastic foam inserts for increased comfort
- Easy-grip tiller adjustment knob
- Rear-view mirror (standard)
- Auto-connecting front-to-rear harness

ADDITIONAL INFORMATION

Weight Capacity	159 kg (350 lbs.)
Maximum Speed	Up to 8.5 km/h (5 mph)
Ground Clearance	8.2 cm (3.2")
Turning Radius	157.5 cm (62")
Overall Length	119 cm (47")
Overall Width	56.5 cm (22")
Front Tyres	9 cm x 25 cm (3.5" x 10") solid
Rear Tyres	10 cm x 25 cm (4" x 10") solid
Suspension	Front and Rear
Range Per Charge	Up to 25 km (15.5 miles)
Total Weight Without Batteries	63.5 kg (140 lbs.)
Heaviest Piece When Disassembled	32 kg (71 lbs.) front section
Standard Seating	Type: Reclining high-back seat with headrest Weight: 13 kg (28 lbs.) Material: Black Vinyl Dimensions: Width: 46 cm (18") Depth: 43 cm (17")
Motor Output	24-volt DC motor, sealed mini-transaxle
Brakes	Electronic, regenerative and electromechanical

(2) 12V deep-cycle

Size: 35AH (standard) or 40 AH (optional)

Weight: 35 AH:15 kg (32.5 lbs.) each

40 AH:15 kg (32.5 lbs.) each

Battery Requirements

Battery Charger

Off-board, 3-amp

Colour

Viper Blue