



**SKU:** S2224

## **PRIDE PATHRIDER 10 DELUXE MOBILITY SCOOTER**

**From: \$3,649.99**

The sleek, sporty Pathrider 10 Deluxe scooter delivers smooth performance and all-new features. With a full complement of advanced standard features including, rugged suspension, a high-visibility automatic front and rear lighting system, feather-touch disassembly, an LED kerb light, and a wraparound easy-drive tiller, the Pathrider 10 Deluxe delivers greater value than ever before.



Scan QR code to  
view product online

## PRODUCT DESCRIPTION

### Key Features

- 159 kg weight capacity
- 157.5 cm turning radius
- Up to 8.5 km/h maximum speed
- Up to 25 km per charge
- Exclusive one-hand feather touch disassembly provides the ultimate quick transport system
- Easy-drive tiller with wraparound handles
- Back-lit battery gauge
- Front and rear suspension
- Lightweight pinch less hinge seat, eight position rotation, sliders, and viscoelastic foam inserts for increased comfort
- Easy-grip tiller adjustment knob
- Rear-view mirror (standard)
- Auto-connecting front-to-rear harness

## ADDITIONAL INFORMATION

<b>Weight Capacity</b>	159 kg (350 lbs.)
<b>Maximum Speed</b>	Up to 8.5 km/h (5 mph)
<b>Ground Clearance</b>	8.2 cm (3.2")
<b>Turning Radius</b>	157.5 cm (62")
<b>Overall Length</b>	119 cm (47")
<b>Overall Width</b>	56.5 cm (22")
<b>Front Tyres</b>	9 cm x 25 cm (3.5" x 10") solid
<b>Rear Tyres</b>	10 cm x 25 cm (4" x 10") solid
<b>Suspension</b>	Front and Rear
<b>Range Per Charge</b>	Up to 25 km (15.5 miles)
<b>Total Weight Without Batteries</b>	63.5 kg (140 lbs.)
<b>Heaviest Piece When Disassembled</b>	32 kg (71 lbs.) front section
<b>Standard Seating</b>	Type: Reclining high-back seat with headrest Weight: 13 kg (28 lbs.) Material: Black Vinyl Dimensions: Width: 46 cm (18") Depth: 43 cm (17")
<b>Motor Output</b>	24-volt DC motor, sealed mini-transaxle
<b>Brakes</b>	Electronic, regenerative and electromechanical

(2) 12V deep-cycle

Size: 35AH (standard) or 40 AH (optional)

Weight: 35 AH:15 kg (32.5 lbs.) each

40 AH:15 kg (32.5 lbs.) each

**Battery Requirements**

**Battery Charger**

Off-board, 3-amp

**Colour**

Viper Blue