

SKU:

PRIDE MOBILITY PATHRIDER 140XL MOBILITY SCOOTER

\$5,999.99

The Pride Pathrider 140XL Mobility Scooter is a premium large scooter with 13" tyres, full suspension and hydraulic brakes. It supports 181 kg and offers smooth, powerful performance for outdoor travel.



Scan QR code to
view product online



VARIATIONS

SKU	Price	Description	Product Colour
S15304	\$5,999.99		Yellow
S2157	\$5,999.99		Red
S15412	\$5,999.99		Silver

PRODUCT DESCRIPTION

The Pathrider 140XL blends strong design with comfort-driven features. It includes full system suspension and sealed hydraulic brakes for better control and smooth handling over rough terrain.

The ergonomic seat offers pillow-top cushioning, swivel base, and adjustable depth, recline, lumbar and headrest settings.

The delta tiller with wraparound handles supports single-hand use.

Dual LED lights, indicators, rear-view mirror and front basket come standard, along with tie-down

points for safe transit.

Who it's for

- For clients wanting a heavy-duty scooter with advanced comfort and braking.
- For therapists seeking a stable, performance-focused scooter with high weight capacity.

Key features

- **Full suspension and hydraulic brakes** – smooth control across uneven terrain
 - **Pillow-top seat with full adjustability** – enhances long-distance comfort
 - **13” pneumatic tyres and 181 kg capacity** – suits various terrains and users
- **Ergonomic tiller and tie-down points** – improves control and secure transport (for the transport of unoccupied scooter)
 - **LED lighting and indicators** – improves visibility day and night

ADDITIONAL INFORMATION

Safe Working Load (SWL)	181 kg
Battery Type	Gel
Range	35 km
Max Speed	10 km/h
Weight	164.6 kg
Suspension	Front & Rear
Product Colour	Red, Yellow, Silver

[Pride Mobility Pathrider 140XL Brochure](#)

[Pride Mobility Pathrider 140XL User Manual](#)

[Pride Mobility Scooter Specification Chart](#)

Talk to us on the phone

Our experienced phone sales team can assist you with any product enquiries and help you make a purchase directly on the phone. We are available anytime between 9am and 5pm, Monday to Friday.

[1300 008 267](tel:1300008267)

WHAT KIND OF BATTERIES DOES THE PATHRIDER 140XL USE?

The batteries are located on the rear section of your scooter. The Pathrider 140XL uses 2 12-volt, deep-cycle 70 Ah batteries that weight 25 kg (55 lbs.) each.

DOES THE PRIDE 140XL PATHRIDER MOBILITY SCOOTER HAVE ANTI TIP WHEELS PLEASE?

140XL Pathrider Mobility Scooter doesn't have anti-tip wheels. 130XL mobility scooter has anti-tip wheels.

IS THERE A WARRANTY WITH PRIDE 140XL PATHRIDER MOBILITY SCOOTER?

Pride 140XL Pathrider Mobility Scooter has two years warranty plus 6 months on batteries.