



**SKU:** S2157

## **PRIDE PATHRIDER 140XL MOBILITY SCOOTER**

**\$5,999.99**

The Pathrider 140 XL Mobility Scooter is an exciting blend of power and precision that features a 24-volt DC motor, full suspension, 33 cm (13") pneumatic tyres and a hydraulic-sealed brake system. Standard luxury touches, including an infinite adjustable wraparound easy drive tiller, high-back pillow-top seat and LED lighting make it the perfect outdoor scooter.



Scan QR code to  
view product online

VARIATIONS

| Image                                                                            | SKU    | Price      | Description | Product Colour  |
|----------------------------------------------------------------------------------|--------|------------|-------------|-----------------|
|  | S15304 | \$5,999.99 |             | Yellow          |
|  | S2157  | \$5,999.99 |             | Candy apple red |
|  | S15412 | \$5,999.99 |             | Silver          |

## PRODUCT DESCRIPTION

### **Pride Pathrider 140XL mobility scooter**

#### **A powerful blend of power and precision for outdoor use**

Designed for outdoor use and on uneven terrain, the Pride Pathrider 140XL is one of our best-selling mobility scooters.

Engineered for powerful performance and exceptional style, key features include full suspension for optimum handling, a hydraulic sealed brake system and large 33cm pneumatic tyres. It is equipped with many standard luxury touches such as the wraparound easy drive tiller and LED lighting.

#### **Key features include:**

- Full suspension for optimum handling
- Sealed hydraulic brakes
- Large 33cm pneumatic tyres for excellent outdoor performance
- Full lighting package including directional signals and headlight
- Easy drive tiller with wraparound handles
- Easily accessible tie down points (for transport of unoccupied scooter)
- Rear view mirror standard
- Front shield standard
- Pillow-top seat

Additional accessories which can be purchased separately:

- Double crutch holder
- Oxygen tank holder
- Rear basket
- Additional rear view mirror

## ADDITIONAL INFORMATION

|                                                                                                            |                                                                                                    |
|------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------|
| <b>Product Safe Working Load (SWL)</b>                                                                     | 181 kg (400 lbs)                                                                                   |
| <b>Product Battery Size / Type</b>                                                                         | 2 x 75 AH Gel                                                                                      |
| <b>Product Range</b>                                                                                       | 35 km                                                                                              |
| <b>Product Speed (Max)</b>                                                                                 | 10 km/hr                                                                                           |
| <b>Product Colour</b>                                                                                      | Candy apple red, Yellow, Silver                                                                    |
| <b>Product Warranty</b>                                                                                    | Mechanical 24 months<br>Electrical 12 months<br>Batteries are covered by a six (6) month warranty. |
| <b>Product Size - in carton(s) (WxLxH) in cm</b>                                                           | 74 x 149 x 46                                                                                      |
| <b>Product Size - in carton(s) (gross weight) in kg</b>                                                    | 125 kg                                                                                             |
| <b>Product Size - overall assembled (WxLxH)* including pulled apart/folded for portable scooters in cm</b> | 69.2 x 142.8 x 106.5                                                                               |
| <b>Product Size - overall assembled (gross weight) in kg</b>                                               | 167.7 kg                                                                                           |
| <b>Product Size - seat base width in cm</b>                                                                | 49.5 cm                                                                                            |
| <b>Product Size - seat base depth in cm</b>                                                                | 44.5 cm                                                                                            |

|                                                         |                                                               |
|---------------------------------------------------------|---------------------------------------------------------------|
| <b>Product Size - seat back height in cm</b>            | 47 cm                                                         |
| <b>Product Size - seat to floor in cm</b>               | 63 cm – 67 cm                                                 |
| <b>Product Size - front/ rear tyre type/size</b>        | 13" / Pneumatic                                               |
| <b>Product Size - front/rear wheel size in inches</b>   | 8"                                                            |
| <b>Product Drive - turning circle in cm</b>             | 210.8 cm                                                      |
| <b>Product Drive - max safe ground clearance in cm</b>  | 13.3 cm                                                       |
| <b>Product Drive - max safe incline slope in degree</b> | 8°                                                            |
| <b>Product Drive - suspension</b>                       | Yes (Front And Rear CTS Suspension With Height Adjustability) |
| <b>Product Charger Type</b>                             | 8 AH                                                          |
| <b>Product Motor Size</b>                               | 4 Pole / 2304 Watt                                            |
| <b>Product Lighting</b>                                 | LED                                                           |
| <b>Product Display</b>                                  | Analog                                                        |