

**SKU:**

# PRIDE MOBILITY PATHRIDER 140XL MOBILITY SCOOTER

**\$5,999.99**

The Pride Pathrider 140XL Mobility Scooter is a premium large scooter with 13" tyres, full suspension and hydraulic brakes. It supports 181 kg and offers smooth, powerful performance for outdoor travel.



Scan QR code to  
view product online

## VARIATIONS

SKU	Price	Description	Product Colour
S15304	\$5,999.99		Yellow
S2157	\$5,999.99		Red
S15412	\$5,999.99		Silver

## PRODUCT DESCRIPTION

The Pathrider 140XL blends strong design with comfort-driven features. It includes full system suspension and sealed hydraulic brakes for better control and smooth handling over rough terrain.

The ergonomic seat offers pillow-top cushioning, swivel base, and adjustable depth, recline, lumbar and headrest settings.

The delta tiller with wraparound handles supports single-hand use.

Dual LED lights, indicators, rear-view mirror and front basket come standard, along with tie-down

points for safe transit.

### Who it's for

- For clients wanting a heavy-duty scooter with advanced comfort and braking.
- For therapists seeking a stable, performance-focused scooter with high weight capacity.

### Key features

- **Full suspension and hydraulic brakes** – smooth control across uneven terrain
  - **Pillow-top seat with full adjustability** – enhances long-distance comfort
  - **13” pneumatic tyres and 181 kg capacity** – suits various terrains and users
- **Ergonomic tiller and tie-down points** – improves control and secure transport (for the transport of unoccupied scooter)
  - **LED lighting and indicators** – improves visibility day and night

## ADDITIONAL INFORMATION

<b>Safe Working Load (SWL)</b>	181 kg
<b>Battery Type</b>	Gel
<b>Range</b>	35 km
<b>Max Speed</b>	10 km/h
<b>Weight</b>	164.6 kg
<b>Suspension</b>	Front & Rear
<b>Product Colour</b>	Red, Yellow, Silver

[Pride Mobility Pathrider 140XL Brochure](#)

[Pride Mobility Pathrider 140XL User Manual](#)

[Pride Mobility Scooter Specification Chart](#)

## **Talk to us on the phone**

Our experienced phone sales team can assist you with any product enquiries and help you make a purchase directly on the phone. We are available anytime between 9am and 5pm, Monday to Friday.

[1300 008 267](tel:1300008267)

### **WHAT KIND OF BATTERIES DOES THE PATHRIDER 140XL USE?**

The batteries are located on the rear section of your scooter. The Pathrider 140XL uses 2 12-volt, deep-cycle 70 Ah batteries that weight 25 kg (55 lbs.) each.

### **DOES THE PRIDE 140XL PATHRIDER MOBILITY SCOOTER HAVE ANTI TIP WHEELS PLEASE?**

140XL Pathrider Mobility Scooter doesn't have anti-tip wheels. 130XL mobility scooter has anti-tip wheels.

### **IS THERE A WARRANTY WITH PRIDE 140XL PATHRIDER MOBILITY SCOOTER?**

Pride 140XL Pathrider Mobility Scooter has two years warranty plus 6 months on batteries.