



**SKU:**

## **PRIDE PATHRIDER 140XL MOBILITY SCOOTER**

**\$5,999.99**

The Pathrider 140 XL Mobility Scooter is an exciting blend of power and precision that features a 24-volt DC motor, full suspension, 33 cm (13") pneumatic tyres and a hydraulic-sealed brake system. Standard luxury touches, including an infinite adjustable wraparound easy drive tiller, high-back pillow-top seat and LED lighting make it the perfect outdoor scooter.



Scan QR code to  
view product online

VARIATIONS

Image	SKU	Price	Description	Product Colour
	S15304	\$5,999.99		Yellow
	S2157	\$5,999.99		Candy apple red
	S15412	\$5,999.99		Silver

## PRODUCT DESCRIPTION

### **Pride Pathrider 140XL mobility scooter**

#### **A powerful blend of power and precision for outdoor use**

Designed for outdoor use and on uneven terrain, the Pride Pathrider 140XL is one of our best-selling mobility scooters.

Engineered for powerful performance and exceptional style, key features include full suspension for optimum handling, a hydraulic sealed brake system and large 33cm pneumatic tyres. It is equipped with many standard luxury touches such as the wraparound easy drive tiller and LED lighting.

#### **Key features include:**

- Full suspension for optimum handling
- Sealed hydraulic brakes
- Large 33cm pneumatic tyres for excellent outdoor performance
- Full lighting package including directional signals and headlight
- Easy drive tiller with wraparound handles
- Easily accessible tie down points (for transport of unoccupied scooter)
- Rear view mirror standard
- Front shield standard
- Pillow-top seat

Additional accessories which can be purchased separately:

- Double crutch holder
- Oxygen tank holder
- Rear basket
- Additional rear view mirror

## ADDITIONAL INFORMATION

<b>Product Safe Working Load (SWL)</b>	181 kg (400 lbs)
<b>Product Battery Size / Type</b>	2 x 75 AH Gel
<b>Product Range</b>	35 km
<b>Product Speed (Max)</b>	10 km/hr
<b>Product Colour</b>	Candy apple red, Yellow, Silver
<b>Product Warranty</b>	Mechanical 24 months Electrical 12 months Batteries are covered by a six (6) month warranty.
<b>Product Size - in carton(s) (WxLxH) in cm</b>	74 x 149 x 46
<b>Product Size - in carton(s) (gross weight) in kg</b>	125 kg
<b>Product Size - overall assembled (WxLxH)* including pulled apart/folded for portable scooters in cm</b>	69.2 x 142.8 x 106.5
<b>Product Size - overall assembled (gross weight) in kg</b>	167.7 kg
<b>Product Size - seat base width in cm</b>	49.5 cm
<b>Product Size - seat base depth in cm</b>	44.5 cm

<b>Product Size - seat back height in cm</b>	47 cm
<b>Product Size - seat to floor in cm</b>	63 cm – 67 cm
<b>Product Size - front/ rear tyre type/size</b>	13" / Pneumatic
<b>Product Size - front/rear wheel size in inches</b>	8"
<b>Product Drive - turning circle in cm</b>	210.8 cm
<b>Product Drive - max safe ground clearance in cm</b>	13.3 cm
<b>Product Drive - max safe incline slope in degree</b>	8°
<b>Product Drive - suspension</b>	Yes (Front And Rear CTS Suspension With Height Adjustability)
<b>Product Charger Type</b>	8 AH
<b>Product Motor Size</b>	4 Pole / 2304 Watt
<b>Product Lighting</b>	LED
<b>Product Display</b>	Analog