



SKU:

PRIDE MOBILITY PATHRIDER 140XL MOBILITY SCOOTER



\$5,999.99

The Pride Pathrider 140XL Mobility Scooter is a premium large scooter with 13" tyres, full suspension and hydraulic brakes. It supports 181 kg and offers smooth, powerful performance for outdoor travel.



Scan QR code to
view product online

VARIATIONS

Image	SKU	Price	Description	Product Colour
	S15304	\$5,999.99		Yellow
	S2157	\$5,999.99		Candy apple red
	S15412	\$5,999.99		Silver

PRODUCT DESCRIPTION

The Pathrider 140XL blends strong design with comfort-driven features. It includes full system suspension and sealed hydraulic brakes for better control and smooth handling over rough terrain. The ergonomic seat offers pillow-top cushioning, swivel base, and adjustable depth, recline, lumbar and headrest settings.

The delta tiller with wraparound handles supports single-hand use.

Dual LED lights, indicators, rear-view mirror and front basket come standard, along with tie-down points for safe transit.

Who it's for

- For clients wanting a heavy-duty scooter with advanced comfort and braking.
- For therapists seeking a stable, performance-focused scooter with high weight capacity.

Key features

- **Full suspension and hydraulic brakes** – smooth control across uneven terrain
- **Pillow-top seat with full adjustability** – enhances long-distance comfort
- **13" pneumatic tyres and 181 kg capacity** – suits various terrains and users
- **Ergonomic tiller and tie-down points** – improves control and secure transport (for the transport of unoccupied scooter)
- **LED lighting and indicators** – improves visibility day and night

ADDITIONAL INFORMATION

Product Colour

Candy apple red, Yellow, Silver