



SKU:

PRIDE MOBILITY PATHRIDER 140XL MOBILITY SCOOTER

\$5,999.99

The Pathrider 140 XL Mobility Scooter is an exciting blend of power and precision that features a 24-volt DC motor, full suspension, 33 cm (13") pneumatic tyres and a hydraulic-sealed brake system. Standard luxury touches, including an infinite adjustable wraparound easy drive tiller, high-back pillow-top seat and LED lighting make it the perfect outdoor scooter.



Scan QR code to
view product online

VARIATIONS

Image	SKU	Price	Description	Product Colour
	S15304	\$5,999.99		Yellow
	S2157	\$5,999.99		Candy apple red
	S15412	\$5,999.99		Silver

PRODUCT DESCRIPTION

Pride Pathrider 140XL mobility scooter

A powerful blend of power and precision for outdoor use

Designed for outdoor use and on uneven terrain, the Pride Pathrider 140XL is one of our best-selling mobility scooters.

Engineered for powerful performance and exceptional style, key features include full suspension for optimum handling, a hydraulic sealed brake system and large 33cm pneumatic tyres. It is equipped with many standard luxury touches such as the wraparound easy drive tiller and LED lighting.

Key features include:

- Full suspension for optimum handling
- Sealed hydraulic brakes
- Large 33cm pneumatic tyres for excellent outdoor performance
- Full lighting package including directional signals and headlight
- Easy drive tiller with wraparound handles
- Easily accessible tie down points (for transport of unoccupied scooter)
- Rear view mirror standard
- Front shield standard
- Pillow-top seat

Additional accessories which can be purchased separately:

- Double crutch holder
- Oxygen tank holder
- Rear basket
- Additional rear view mirror

ADDITIONAL INFORMATION

Product Safe Working Load (SWL)	181 kg (400 lbs)
Product Battery Size / Type	2 x 75 AH Gel
Product Range	35 km
Product Speed (Max)	10 km/hr
Product Colour	Candy apple red, Yellow, Silver
Product Warranty	Mechanical 24 months Electrical 12 months Batteries are covered by a six (6) month warranty.
Product Size - in carton(s) (WxLxH) in cm	74 x 149 x 46
Product Size - in carton(s) (gross weight) in kg	125 kg
Product Size - overall assembled (WxLxH)* including pulled apart/folded for portable scooters in cm	69.2 x 142.8 x 106.5
Product Size - overall assembled (gross weight) in kg	167.7 kg
Product Size - seat base width in cm	49.5 cm
Product Size - seat base depth in cm	44.5 cm

Product Size - seat back height in cm	47 cm
Product Size - seat to floor in cm	63 cm – 67 cm
Product Size - front/ rear tyre type/size	13" / Pneumatic
Product Size - front/rear wheel size in inches	8"
Product Drive - turning circle in cm	210.8 cm
Product Drive - max safe ground clearance in cm	13.3 cm
Product Drive - max safe incline slope in degree	8°
Product Drive - suspension	Yes (Front And Rear CTS Suspension With Height Adjustability)
Product Charger Type	8 AH
Product Motor Size	4 Pole / 2304 Watt
Product Lighting	LED
Product Display	Analog