


SKU: S10833

PRIDE MOBILITY S19 QUEST MOBILITY SCOOTER

\$3,899.99

The Pride S19 Quest Mobility Scooter is a lightweight folding scooter ideal for travel. With a lithium battery, LED lighting and a 19 km range, it supports both indoor and outdoor mobility.



Scan QR code to
view product online



PRODUCT DESCRIPTION

The Pride S19 Quest Mobility Scooter is a practical folding scooter designed for those needing a compact, transportable mobility scooter. It folds down in three easy steps into a suitcase-like shape and can be stored upright or flat – ideal for small spaces or car boots.

Powered by a lithium-ion battery, it travels up to 19 km at speeds of up to 6 km/h, supporting daily activities with ease.

It has an ultra-light design, with an adjustable tiller to allow for tailored control. Two extra assist wheels improve stability and handling.

Front and rear LED lights enhance visibility, making the S19 suitable for both indoor and outdoor use. Its compact footprint makes it easy to manoeuvre in tight spaces, without compromising on safety or convenience.

Who it's for

- For clients looking for a highly portable, folding scooter for everyday and travel use.
- For therapists recommending a lightweight, low-maintenance scooter with stable handling.

Key features

- **Folds in 3 simple steps** – enables easy transport and compact storage
- **Lithium-ion battery with 19 km range** – supports daily travel needs
- **Lightweight, suitcase-style foldable frame** – ideal for travel and small storage spaces
- **Stows upright or laid flat** – suits versatile storage needs
- **LED front and rear lighting** – improves safety in low-light conditions
- **Adjustable tiller and assist wheels** – support handling and personal control

ADDITIONAL INFORMATION

Product Safe Working Load (SWL)	115 kg (253.5 lbs)
Product Battery Size / Type	11.5 AH Lithium
Product Range	20 km
Product Speed (Max)	6 km/h
Product Colour	White
Product Warranty	Mechanical 12 months Electrical 12 months Batteries are covered by a six (6) month warranty.
Product Size - in carton(s) (WxLxH) in cm	46 x 51 x 88
Product Size - in carton(s) (gross weight) in kg	34 kg
Product Size - overall assembled (WxLxH)* including pulled apart/folded for portable scooters in cm	49.5 x 92.7 x 96
Product Size - overall assembled (gross weight) in kg	25.8 kg
Product Size - seat base width in cm	41 cm
Product Size - seat base depth in cm	33 cm

Product Size - seat back width in cm	41 cm
Product Size - seat back height in cm	30 cm
Product Size - seat to floor in cm	51 cm
Product Size - front/ rear tyre type/size	7" / Solid
Product Size - front/rear wheel size in inches	5"
Product Drive - turning circle in cm	164 cm
Product Drive - max safe ground clearance in cm	6 cm
Product Drive - max safe incline slope in degree	6°
Product Drive - suspension	No
Product Charger Type	2.5 AH
Product Motor Size	2 Pole / 864 Watt
Product Lighting	LED
Product Display	Analog