

SKU: S12062



## PROTENS & EMS DUAL UNIT WITH PREMIUM DISPLAY

~~\$199.99~~ **\$139.99**

Digital Tens and EMS With Hard Case. Combined EMS and TENS Unit.

Digital Tens has the ability to store your treatment times and schedules which can be examined by your chiropractor. Unit comes with TENS stimulator unit, carrying case, 2X TENS leads, 3XAAA batteries, operating manual and packet of 4 Electrodes. Includes a belt clip.

**5 Years Unlimited Warranty**

**You may be entitled to receive a 100% rebate from your health fund - Subject to your health cover**



Scan QR code to  
view product online

## **PRODUCT DESCRIPTION**

ILS Protens & EMS Dual Unit has fully adjustable parameters and provides relief from acute or chronic pain such as back pain, sports injury, general muscle pain, etc. Battery powered electrical unit which uses electrodes to deliver electrical pulses to the nerve fibres relieving pain.

## ADDITIONAL INFORMATION

|                          |                                                       |
|--------------------------|-------------------------------------------------------|
| <b>Brand</b>             | ILS physio                                            |
| <b>System Type</b>       | TENS and EMS                                          |
| <b>Level</b>             | 0 - 25 Level                                          |
| <b>Modes</b>             | TENS□Program 8<br>EMS□Program 6                       |
| <b>Pulse Intensity</b>   | Adjustable 0-100mA peak into 500ohm load each channel |
| <b>Pulse width range</b> | TENS□156~260 uS (pk)<br>EMS□200~300 (pk)              |
| <b>Intensity</b>         | 0~80mA                                                |
| <b>Ramp Time</b>         | 1 - 30 seconds                                        |
| <b>Power indicator</b>   | LCD                                                   |
| <b>Cycle Off Time</b>    | 1 - 60 seconds                                        |
| <b>Timer</b>             | 5~60 min (5 minutes for one unit)                     |
| <b>Power Source</b>      | 4.5V Battery (AAA×3)                                  |
| <b>Channel Type</b>      | Dual                                                  |