

SKU: S12062



PROTENS & EMS DUAL UNIT WITH PREMIUM DISPLAY

~~\$199.99~~ **\$139.99**

Digital Tens and EMS With Hard Case. Combined EMS and TENS Unit.

Digital Tens has the ability to store your treatment times and schedules which can be examined by your chiropractor. Unit comes with TENS stimulator unit, carrying case, 2X TENS leads, 3XAAA batteries, operating manual and packet of 4 Electrodes. Includes a belt clip.

5 Years Unlimited Warranty

You may be entitled to receive a 100% rebate from your health fund - Subject to your health cover



Scan QR code to
view product online

PRODUCT DESCRIPTION

ILS Protens & EMS Dual Unit has fully adjustable parameters and provides relief from acute or chronic pain such as back pain, sports injury, general muscle pain, etc. Battery powered electrical unit which uses electrodes to deliver electrical pulses to the nerve fibres relieving pain.

ADDITIONAL INFORMATION

Brand	ILS physio
System Type	TENS and EMS
Level	0 - 25 Level
Modes	TENS□Program 8 EMS□Program 6
Pulse Intensity	Adjustable 0-100mA peak into 500ohm load each channel
Pulse width range	TENS□156~260 uS (pk) EMS□200~300 (pk)
Intensity	0~80mA
Ramp Time	1 - 30 seconds
Power indicator	LCD
Cycle Off Time	1 - 60 seconds
Timer	5~60 min (5 minutes for one unit)
Power Source	4.5V Battery (AAA×3)
Channel Type	Dual