




SKU:

Q700 R SEDEO ERGO POWER CHAIR

Discover the ultimate driving experience with the striking Q700® R, the high-performance rear-wheel drive power wheelchair. Innovative Spirit Level technology and all-round suspension ensures optimal comfort on all four wheels. It can climb kerbs of up to 100 mm with a kerb climber or 80 mm without it. Enjoy comfort as standard with the amazing SEDEO® ERGO, the biomechanical seating system that mirrors your body's natural position for an ideal seated position at all times, even during transitions! Plus, stay connected with the world around you with a selection of easy-to-use controls, personalise your drive with QUICKIE's patented assignable buttons and 6 seated memory positions.



Scan QR code to
view product online

PRODUCT DESCRIPTION

SEDEO ERGO - INTELLIGENT SEATING

It's the first-class seat that's designed specifically to fit you. The SEDEO ERGO with patented Biometric repositioning is more effective than any other standing wheelchair seating system when it comes to independent repositioning throughout the day.

REAR-WHEEL POWERED WHEELCHAIR

Enjoy the ultimate driving experience with the Q700 R SEDEO ERGO, QUICKIE's most advanced rear-wheel drive powered wheelchair! Featuring Spirit Level technology, the advanced SEDEO ERGO biomechanical seating system and innovative controls.

ADDITIONAL INFORMATION

Seat Width:	400 – 560 mm
Seat Height:	430 – 480 mm
Seat Depth:	400 – 560 mm
Backrest Height:	510 – 610 mm
Overall Width:	640 mm
Overall Length:	Max. 1180-1370 mm
Speed:	6 kph high torque, 10 kph standard
Battery Size:	60 Ah, 72 Ah optional
Max. Range:	40 km
Seat Tilt:	30° / 50°
Turning Diameter:	950 mm
Max. Safe Slope:	10°
Max Kerb Climb:	100 mm with Kerb Climber / 80 mm without KC
Electronics:	R-Net controller; lights & indicators
Max. User Weight:	160 kg
Brand	Sunrise Medical