


SKU:

HERO MEDICAL QUAD SEAT WALKER - 8 INCH WHEELS

\$199.99



The 8 inch Medical quad Seat Walker is made out of heavy duty construction with a safe working limit of 130kg and is Australian standards certified. The new “light touch” lock brake system and 20cm flat ridged wheels minimise wheel wobble and increase grip. The easy fold rounded backrest offers comfort, and is ideal for indoor and outdoor use. This seat walker can be folded without removing the storage bag (which is made out of nylon and prevents visibility of valuable items. e.g. wallet, handbag, keys).



Scan QR code to
view product online

VARIATIONS

Image	SKU	Price	Description	Colour
 A blue four-wheeled walker with a black seat and handgrips.	S11173	\$199.99		Blue
 A burgundy four-wheeled walker with a black seat and handgrips.	S11172	\$199.99		Burgundy

PRODUCT DESCRIPTION

The Freedom Module 8 inch Seat Walker comes with the following features:

- Ergonomic hand grips
- Height adjustable handles
- Lightweight aluminium construction
- Light-touch lock brakes
- Easy fold rounded padded back rest
- Spill proof vinyl padded seat
- Nylon bag prevents visibility of valuable items
- 20cm flat ridged wheels

ADDITIONAL INFORMATION

Brand	Hero Medical
Weight Capacity / Max SWL	140kg (308 lbs)
Overall height (floor to handles)	Min 83cm - Max 94cm
Overall length	69cm
Overall width	62cm
Seat depth	33cm
Seating width	44cm
Seat height	56cm
Folded width	23cm
Outer width (Including back castors)	58cm
Unit weight	8.4kg
Castor size	20 cm / 8 in
Warranty	12 Months Warranty. Please note: Brake maintenance and adjustment is not covered under warranty