

SKU:

## **HERO MEDICAL QUAD SEAT WALKER - 8 INCH WHEELS**

**\$199.99**



The 8 inch Medical quad Seat Walker is made out of heavy duty construction with a safe working limit of 130kg and is Australian standards certified. The new “light touch” lock brake system and 20cm flat ridged wheels minimise wheel wobble and increase grip. The easy fold rounded backrest offers comfort, and is ideal for indoor and outdoor use. This seat walker can be folded without removing the storage bag (which is made out of nylon and prevents visibility of valuable items. e.g. wallet, handbag, keys).



Scan QR code to  
view product online

### VARIATIONS

Image	SKU	Price	Description	Colour
 A four-wheeled walker with a blue frame, black seat, and black handgrips.	S11173	\$199.99		Blue
 A four-wheeled walker with a burgundy frame, black seat, and black handgrips.	S11172	\$199.99		Burgundy

## **PRODUCT DESCRIPTION**

The Freedom Module 8 inch Seat Walker comes with the following features:

- Ergonomic hand grips
- Height adjustable handles
- Lightweight aluminium construction
- Light-touch lock brakes
- Easy fold rounded padded back rest
- Spill proof vinyl padded seat
- Nylon bag prevents visibility of valuable items
- 20cm flat ridged wheels

## ADDITIONAL INFORMATION

<b>Brand</b>	Hero Medical
<b>Weight Capacity / Max SWL</b>	140kg (308 lbs)
<b>Overall height (floor to handles)</b>	Min 83cm - Max 94cm
<b>Overall length</b>	69cm
<b>Overall width</b>	62cm
<b>Seat depth</b>	33cm
<b>Seating width</b>	44cm
<b>Seat height</b>	56cm
<b>Folded width</b>	23cm
<b>Outer width (Including back castors)</b>	58cm
<b>Unit weight</b>	8.4kg
<b>Castor size</b>	20 cm / 8 in
<b>Warranty</b>	12 Months Warranty. Please note: Brake maintenance and adjustment is not covered under warranty