


SKU:

HERO MEDICAL 8" STANDARD SEAT WALKER

\$199.99



The Hero Medical 8" Standard Seat Walker supports everyday walking with added comfort and safety. Featuring large 8" castors, a folding frame and cushioned seating, it's suitable for both indoor and outdoor use.



Scan QR code to
view product online

VARIATIONS

SKU	Price	Description	Colours
S11173	\$199.99	The Hero Medical 8" Standard Seat Walker supports everyday walking with added comfort and safety. Featuring large 8" castors, a folding frame and cushioned seating, it's suitable for both indoor and outdoor use.	Blue
S11172	\$199.99	The Hero Medical 8" Standard Seat Walker supports everyday walking with added comfort and safety. Featuring large 8" castors, a folding frame and cushioned seating, it's suitable for both indoor and outdoor use.	Burgundy

GALLERY IMAGES



PRODUCT DESCRIPTION

The Hero Medical 8" Standard Seat Walker is a practical mobility aid for individuals who need consistent walking support with the ability to rest as needed. Its 8" flat-ridged castors provide improved stability and smoother travel across varied indoor and outdoor surfaces.

The walker includes a padded seat and backrest for rest stops, ergonomic height-adjustable handles for optimal walking posture, and light-touch loop brakes with a locking function for added control.

Its lightweight aluminium frame folds easily for transport and storage. A built-in under-seat bag keeps essentials close during daily activities.

Who it's for

- For clients who require a cost-effective walking aid with larger castors for smoother indoor and outdoor mobility.
- For therapists prescribing a compact, adjustable walker that promotes safe posture and supports rest.

Key features

- **8" castors** – provides smooth manoeuvrability across indoor and outdoor environments
 - **Lightweight aluminium frame** – easy to transport and store
 - **Four-wheeled design** – supports stable, everyday walking
 - **Padded seat and backrest** – offers comfort during breaks
 - **Height-adjustable handles** – supports ergonomic positioning
 - **Light-touch loop brakes** – improves safety and control with locking option
- **Folding frame with under-seat storage bag** – keeps personal items secure and accessible

ADDITIONAL INFORMATION

Safe working load (SWL)	140 kg
Gross Weight	8.4 kg
Seat Height	56 cm
Seat Width	44 cm
Handle Height	83 cm – 94 cm
Outer Width	58 cm
Front castor	8"
Rear castor	8"
Frame material	Aluminium
Brand	Hero Medical

[Hero Medical 8" Standard Seat Walker User Manual](#)

Talk to us on the phone

Our experienced phone sales team can assist you with any product enquiries and help you make a purchase directly on the phone. We are available anytime between 9am and 5pm, Monday to Friday.

[1300 008 267](tel:1300008267)