



#### SKU:

# HERO MEDICAL 8" STANDARD SEAT WALKER

\$199.99

The Hero Medical 8" Standard Seat Walker supports everyday walking with added comfort and safety. Featuring large 8" castors, a folding frame and cushioned seating, it's suitable for both indoor and outdoor use.



Scan QR code to view product online



## **VARIATIONS**

Image	SKU	Price	Description	Colours
×	S11173	\$199.99		Blue
×	S11172	\$199.99		Burgundy



#### PRODUCT DESCRIPTION

The Hero Medical 8" Standard Seat Walker is a practical mobility aid for individuals who need consistent walking support with the ability to rest as needed. Its 8" flat-ridged castors provide improved stability and smoother travel across varied indoor and outdoor surfaces.

The walker includes a padded seat and backrest for rest stops, ergonomic height-adjustable handles for optimal walking posture, and light-touch loop brakes with a locking function for added control. Its lightweight aluminium frame folds easily for transport and storage. A built-in under-seat bag keeps essentials close during daily activities.

#### Who it's for

- For clients who require a cost-effective walking aid with larger castors for smoother indoor and outdoor mobility.
- For therapists prescribing a compact, adjustable walker that promotes safe posture and supports rest.

#### **Key features**

- 8" castors provides smooth manoeuvrability across indoor and outdoor environments
- Lightweight aluminium frame easy to transport and store
- Four-wheeled design supports stable, everyday walking
- Padded seat and backrest offers comfort during breaks
- Height-adjustable handles supports ergonomic positioning
- Light-touch loop brakes improves safety and control with locking option
- Folding frame with under-seat storage bag keeps personal items secure and accessible



### **ADDITIONAL INFORMATION**

Safe working load (SWL) 140 kg

**Gross Weight** 8.4 kg

**Seat Height** 56 cm

Seat Width 44 cm

**Handle Height** 83 cm - 94 cm

Outer Width 58 cm

Front castor 8"

**Rear castor** 8"

Frame material Aluminium

Brand Hero Medical