



**SKU:**

## **WALKING STICK QUAD CANE WITH SMALL BASE**

**\$69.99**

The small base quad cane features an adjustable height to suit individual user requirements. The swan neck handle distributes weight evenly while the strong handgrip enhances user comfort.



Scan QR code to  
view product online

**VARIATIONS**

Image	SKU	Price	Description	Colour
	S062	\$69.99		Bronze
	S9041	\$69.99		Silver

## **PRODUCT DESCRIPTION**

The quad cane's handle rotates for left or right hand use and the low profile offset base does not obstruct user's step. A locking collar prevents rattles and a combination lightweight aluminium and steel construction decreases weight while at the same time increasing strength.

## **ADDITIONAL INFORMATION**

<b>Tip Size</b>	13 mm
<b>Base Width</b>	140 mm
<b>Base Depth</b>	210 mm
<b>Height Adjustable</b>	750 mm - 980 mm
<b>Max SWL</b>	100 kg
<b>Brand</b>	Hero Medical