




SKU:

QUANTUM 4FRONT2 HD

The 4Front® is a quiet, more responsive front-wheel drive power chair that features automotive-grade suspension with unprecedented comfort and ride quality. The 4Front's CASE (Caster Angle Sensor Encoder) technology and SRS (Smooth Ride Suspension) provide stability for smooth obstacle transitions and unmatched driving performance. The 4Front features standard speeds of up to 6 mph in the lowered position, and speeds up to 3.2 mph while elevated at 10" with our optional safe seat elevation. Standard LED fender lights allow you to see and be seen.



Scan QR code to
view product online

PRODUCT DESCRIPTION

- CASE (Caster Angle Sensor Encoder) technology
- Unmatched driving performance and ride quality
- Advanced stability for smooth obstacle transitions
- Superior, automotive-grade SRS (Smooth Ride Suspension)
- Q-Logic 3 EX Advanced Drive Controls
- LED front and rear fender lights
- Quiet ride technology
- Available with WC19 crash-test approved, occupied transit option with all configurations

ADDITIONAL INFORMATION

Brand	Quantum
Type	Heavy Duty
Drive Wheel Configuration	Front Wheel Drive
Weight Capacity	204 kg
Transit Approved	Yes
Max Speed	8 km/h
Base Width	24"
Base Length	40"
Pediatric Seating	No
Seat Width Range	20" - 24"
Seat Depth Range	20" - 24"
Seat To Floor Height Range	18" - 20"
Posterior Tilt Option	0° - 50°
Anterior Tilt Option	No
Recline Option	85° - 175°
Legrest Elevation Option	Yes
Seat Elevation Option	10"
Power Standing Option	No