




QUANTUM Q6 EDGE 3 POWER CHAIR

With features like standard 4-pole motors, Mid-Wheel 6® drive design and ATX Suspension, the Q6 Edge® is engineered to meet the performance needs of the most active user. The Q6 Edge accepts a complete range of seating and electronics options to deliver exceptional rehab capability at an excellent value.

Upgraded SRS (Smooth Ride Suspension)
Optional 4.5 mph with iLevel®
Optional iLevel's patented Extra Stability Technology® enhances safety while performing transfers, reach, and numerous other activities of daily living
iLevel is retrofittable on the Edge 3
Fender lights, USB charger, and drive wheel color accents included
7 bold new colors and a total of 19 color options
Made in America of foreign and domestic parts
Available with WC19 crash-test approved occupied transit kit option with captain's seating, Synergy® seating, and TRU-Balance® 3 seating configurations.



SKU: N/A

Categories: [Adult](#), [Mid Wheel Drive](#), [Power Wheelchairs](#), [Powered Chair](#), [Quantum](#)

GALLERY IMAGES



PRODUCT DESCRIPTION

The design of the Edge 3 is for the most advanced power chair experience ever! Upgraded SRS (Smooth Ride Suspension) provides a smooth, comfortable ride and advanced stability. Optional iLevel® is available at 4.5 mph with up to 12" of power adjustable seat height. Standard features include front and rear LED fender lights, a USB charger, and drive wheel color accents. The Edge 3 is

available in a total of 19 color options, including seven brand-new, bold colors.

ADDITIONAL INFORMATION

Drive Wheels	Standard: 14" Solid Optional: 14" Pneumatic
Caster Wheels	Front: 6" Solid Rear: 6" Solid
Maximum Speed	Up to 5 mph (Optional: 6 mph)
Ground Clearance	3" (frame center)
Turning Radius	20"
Overall Length	35.5" without foot riggings
Base Width	24"
Tilt/Recline Manual	Yes
Battery Size	NF-22
Battery Weight	38 lbs
Available Electronics	75A Q-Logic 2 NE 75A Q-Logic 2 NE+ 100A Q-Logic 2 EX
Battery Charger	8A Off-board
Motor Packages	4-Pole 5 mph ¹ or 6 mph ²
Weight Capacity	300 lbs
Base Weight	138 lbs
Brand	Quantum