




SKU:

QUANTUM EDGE HD

The Quantum Edge HD power wheelchair is designed for users who need a high weight capacity and excellent drive quality. With a weight capacity of up to 204 kg and a superior suspension system, this chair provides enhanced traction and stability.



Scan QR code to
view product online

PRODUCT DESCRIPTION

The Quantum Edge HD power wheelchair is an excellent choice for users who require a high weight capacity and advanced drive quality. With a weight capacity of up to 204 kg, the chair is built to accommodate larger individuals while still providing a smooth and stable ride. The Edge HD features a superior suspension system that delivers enhanced traction performance and stability, making it an ideal choice for both indoor and outdoor use. In addition to its weight capacity and suspension system, the Edge HD offers the option of iLevel technology, providing speed and height for a superior driving experience. Overall, a top-of-the-line power wheelchair.

ADDITIONAL INFORMATION

Brand	Quantum
Type	Heavy Duty
Drive Wheel Configuration	Mid Wheel Drive
Weight Capacity	204 kg
Transit Approved	Yes
Max Speed	8 km/h
Base Width	26.5"
Base Length	36"
Pediatric Seating	No
Seat Width Range	16" - 28"
Seat Depth Range	16" - 28"
Seat To Floor Height Range	18.5" - 20.5"
Posterior Tilt Option	0° - 50°
Anterior Tilt Option	No
Recline Option	85° - 175°
Legrest Elevation Option	Yes
Seat Elevation Option	10"
Power Standing Option	No