



SKU:

QUANTUM R44 POWER CHAIR

Striking Performance Rear-wheel Drive

The new R44 rear-wheel drive power chair features a high-forward pivot long trail arm, a state-of-the-art design that provides dramatically improved obstacle climbing performance and a smoother ride. The R44 also has outstanding tight space manoeuvrability, a sliding battery tray for easy access and accepts rehab seating and electronics. With its stylish, sleek design and a standard speed of up to 9.65 km/h (6mph), the R44 is a superior power chair for even the most active clients.



Scan QR code to
view product online

PRODUCT DESCRIPTION

- Narrow 64.1 cm (25") width and compact 88 cm length
- Outstanding tight space manoeuvrability
- Maximum speed up to 9.65 km/h (6 mph)
- High-performance, high-efficiency motor package
- Accepts rehab seating
- Swing-away joystick
- Accommodates full range of Pride® seating and accessories
- Sliding battery tray for easy access

ADDITIONAL INFORMATION

DRIVE WHEELS	35.5 cm (14") (pneumatic or solid available)
CASTERS	Front: 23 cm (9") or 25 cm (10")
ANTI-TIPS	6.4 cm (2.5") solid (rear)
SUSPENSION	Active-Trac® ATX Suspension (with long trail arm)
MAXIMUM SPEED	Up to 9.65 km/h (6 mph)
BRAKES	Intelligent braking (electronic, regenerative disc brakes)
GROUND CLEARANCE	8.25 cm (3.25") at battery tray
TURNING RADIUS	65.4 cm (25.75")
OVERALL SIZE	Length: 110 cm (43.25") 5/ 88 cm (35") 6 Width: 64.1 cm (25") Synergy® Available widths: 25 cm – 56 cm (10"-22") Available depths: 25 cm – 51 cm (10"-20") TRU-Balance® Power Positioning System Adjustable Width: 35.5 cm – 56 cm (14"-22") Adjustable Depth: 35.5 cm – 56 cm (14"-20")
SEATING	Euro 2 Available widths: 43 cm – 51 cm (17"-20") Available depths: 43 cm – 51 cm (17"-20") High-back Seat (22" x 20" max. dim.) Medium-back Seat (22" x 20" max. dim.)
DRIVETRAIN	Two-motor, rear-wheel drive
BATTERIES	Two 12 volt, deep cycle 70 AH (Group 24)
RANGE PER CHARGE	Variable up to 31.38 km (19.5 miles)

BATTERY CHARGER	8A, off-board
WEIGHT CAPACITY	136 kg (300 lbs.) maximum (user weight)
BASE WEIGHT	58.51 kg (129 lbs.)
BATTERY WEIGHT	23.13 kg (51 lbs.) (each)
BRAND	Quantum