



**SKU:** S15784

## **QUICKIE Q100 POWER CHAIR**

**\$4,995.99**

**FREE BONUS HEADREST Valued at  
\$499.99**

### **THE ULTRA-COMPACT INDOOR / OUTDOOR POWERCHAIR.**

From SMART base development the Q100 R powerchair boasts an ultra-small footprint of 540 mm wide and 1010 mm long, with a tiny 1600 mm turning circle, so navigating even the tightest of spaces is a breeze.

Going out? The QUICKIE Q100 R provides a stable ride on urban roads and paths, ideal for those day trips out.

Unbelievably light at 95 kg, the Q100 R is also easy to hoist in and out of your car.

### **THE SMALL POWERCHAIR THAT'S BIG ON GOING OUT**

Popping out? The Q100 R's indoor capabilities are just as at home outdoors. Painstakingly researched, its compact dimensions provide highly-deceptive traction and stability over a variety of terrains, hills and even when tackling kerbs.

### **SIT BACK AND RELAX**

The Q100 R's seating has been designed to fit you. With plenty of settings and a 5-strap

tension adjustable back to meet postural requirements, take comfort in knowing that if your needs change, so can your seat.

Seat adjustments include:

Seat width from 430-480 mm reduced to 380 mm  
with an optional in-rig kit  
Seat depth from 405 - 510 mm  
Optional Legrest width from 15 - 19"  
Backrest angle from 0, 4, 8 and 12 degrees  
Seat angle from 0, 3 and 6 degrees  
Side guard from 0 - 50 mm  
Optional seat height from 450 - 490 mm  
Armrest height 230 - 305 mm



Scan QR code to  
view product online

## **PRODUCT DESCRIPTION**

QUICKIE Q100 power chair uses an innovative testing method called SMART base technology. This helped us to understand how geometry changes to a drive base can affect performance, creating a precision-measured TRUE indoor/outdoor powerchair.

## ADDITIONAL INFORMATION

|                               |                                    |
|-------------------------------|------------------------------------|
| <b>Brand</b>                  | Sunrise Medical                    |
| <b>Seat Width</b>             | 430 – 500 mm (17" – 19")           |
| <b>Seat Height</b>            | 450 mm (17.5")                     |
| <b>Seat Depth</b>             | 405 – 505 mm (16" – 20")           |
| <b>Backrest Height</b>        | 500 mm (20")                       |
| <b>Overall Width</b>          | 540 mm (21")                       |
| <b>Overall Length:</b>        | 1010 mm (39.5")                    |
| <b>Speed</b>                  | 6 kph                              |
| <b>Battery Size</b>           | 55 Ah                              |
| <b>Max. Range</b>             | 31 km (19 miles)                   |
| <b>Seat Tilt</b>              | 0° – 6°                            |
| <b>Back Recline</b>           | 0°, 4°, 8° & 12° (manual)          |
| <b>Turning Diameter</b>       | 1600 mm (63")                      |
| <b>Max. Safe Slope</b>        | 6°                                 |
| <b>Max Kerb Climb</b>         | 70 mm – 100 mm (with Kerb Climber) |
| <b>Electronics</b>            | VR2                                |
| <b>Min. Wheelchair Weight</b> | from 95 kg (incl. 55 Ah batteries) |
| <b>Max. User Weight</b>       | 125 kg                             |
| <b>Colour</b>                 | Blue                               |