



SUNRISE MEDICAL QXI SCRIPTED MANUAL CHAIR

The Quickie Advantage Family includes QXi, a customisable lightweight folding chair and QRi, a rigid wheelchair designed for the long term, semi active or active user.

They offer plenty of adjustability and options, including a transport package, to allow these chairs to be tailored to meet a wide range of common user needs.



SKU: N/A

Category: [Wheelchairs - Largest Wheelchair Selection for folding, Transit and Manual](#)

GALLERY IMAGES



PRODUCT DESCRIPTION

Frame Style	Folding - Traditional
Transport Weight	8 kg
Overall Weight	13 kg
User Weight Capacity	120 kg
Seat Width	14" to 20"
Seat Depth	14" to 20"
Back Angle Adjustment	-15° to 3°
Front Seat to Floor Height	12.5" to 19.5"
Rear Seat to Floor Height	12" to 19"
Hanger Angle & Elevating	60°, 70°, 80°
Warranty	Limited Lifetime

ADDITIONAL INFORMATION

Brand Sunrise Medical