



SKU:


## INVACARE REA AZALEA TILT IN SPACE WHEELCHAIR

~~\$4,999.99~~ \$4,999.99

The Rea Azalea is the natural choice for the moderately active user as well as for the more passive one.

The comfortable seating system allows the user to alter the seat position, always ensuring maximum comfort.

Combined with good driving characteristics makes the Azalea the ideal choice. The Azalea is available in either Self-Propelling or Transit options.



Scan QR code to  
view product online

**VARIATIONS**








Image	SKU	Price	Description	Wheelchair Type	Size
	S10366	\$4,999.99		Transit	14"
	S9688	\$4,699.99		Self Propelled	16"
	S9687	\$4,699.99		Transit	16"
	S9524	\$4,699.99		Self Propelled	18"

Image	SKU	Price	Description	Wheelchair Type	Size
	S9414	\$4,699.99		Transit	18"
	S9736	\$4,699.99		Self Propelled	20"
	S9735	\$4,699.99		Transit	20"

## PRODUCT DESCRIPTION

The Invacare Rea Azalea is a quality and durable Tilt n' Space Wheelchair option, available in varying seat widths in either a self-propelled or transit version. Comfortable padded seating, backrest, armrests and headrest ensures optimal seating comfort.

The Azalea offers a 25° seat tilt, and 30° backrest recline, adjusted with ease using a manual lever. Backrest has a biomechanical positioning hinge which follows the user's natural back movement when adjusting the backrest angle. The backrest plate is a strong and sturdy plate, offering the client stability whether upright or in tilt/recline. Multi adjustable backrest can be adjusted in height and in width for different body shapes.

### Key Features

**Tilt-in-Space** - The Rea Azalea allows users to easily adjust their seating position with its tilt-in-space feature. This helps to reduce pressure sores, improve posture, and provide comfort by redistributing weight and supporting the user's body more effectively.

**Customisable Seating** - This wheelchair offers multiple seating options, including adjustable backrests, seat depths, seat widths, headrest, armrest and legrest. The system allows for personalised configuration to match the unique needs of the user, ensuring maximum comfort and optimal positioning.

**Durable Design** - Built with durability and stability in mind, the Rea Azalea's frame is strong yet lightweight. Its robust design supports weights up to 135 kg and ensures long-term performance even under continuous use.

**Smooth Manoeuvrability** - The chair features high-quality wheels and casters that allow for smooth and easy movement, whether indoors or outdoors. This ensures that caregivers can assist with transfers and positioning without difficulty. Turning radius of 77 cm.

**Easy Adjustments and Maintenance** - The Rea Azalea is designed with ease of use in mind. The intuitive tilt mechanism, tool-free adjustments, and straightforward maintenance make it easy for caregivers to provide necessary care and adjustments without complications. Easy to wipe down and clean multi-stretch PU coated fabric.

## ADDITIONAL INFORMATION

<b>Brand</b>	Invacare
<b>Max SWL</b>	135 kg
<b>Seat depth</b>	43cm – 50cm
<b>Seat width</b>	37cm – 59cm depending on model
<b>Model</b>	14" Transit, 16" Self Propelled, 16" Transit, 18" Self Propelled, 18" Transit, 20" Self Propelled, 20" Transit
<b>Colour</b>	Pearl Grey