



SKU:

## INVACARE REA AZALEA TILT IN SPACE WHEELCHAIR

~~\$4,999.99~~ \$4,999.99

The Invacare Azalea Tilt-in-Space Wheelchair supports individuals with complex seating needs through customisable tilt and recline functions. With multiple adjustments and pressure care features, it delivers tailored comfort and clinical stability.



Scan QR code to  
view product online

**VARIATIONS**








Image	SKU	Price	Description	Wheelchair Type	Size
	S10366	\$4,999.99		Transit	14"
	S9688	\$4,699.99		Self Propelled	16"
	S9687	\$4,699.99		Transit	16"
	S9524	\$4,699.99		Self Propelled	18"

Image	SKU	Price	Description	Wheelchair Type	Size
	S9414	\$4,699.99		Transit	18"
	S9736	\$4,699.99		Self Propelled	20"
	S9735	\$4,699.99		Transit	20"

## PRODUCT DESCRIPTION

The Invacare Azalea Tilt-in-Space Wheelchair is designed to meet postural, pressure care and mobility needs for people with higher support requirements.

Available in self-propelled and transit versions and various seating widths, the Azalea caters to a wide range of body shapes and postural needs.

The Azalea provides up to 25° of seat tilt and 30° of backrest recline, both adjusted via a manual lever. A biomechanical backrest hinge moves with the person's body for more natural support during repositioning. The Dual Stability System (DSS) shifts the seat system inside the frame during tilt to reduce shear and maintain balance.

Adjustability is a core feature: seat depth, backrest angle, armrest height and legrest length can all be tailored. The system also includes an adjustable headrest, lateral supports, calf pads and a seatbelt as standard.

The durable frame supports up to 135 kg, and the wheelchair remains easy to manoeuvre indoors or out, with a 77 cm turning radius.

Its multi-stretch PU fabric is easy to clean, while tool-free adjustments simplify ongoing setup and care.

### Who it's for

- For clients requiring tilt-in-space functionality for postural stability and pressure redistribution.
- For therapists prescribing adjustable, high-support manual mobility for clients with complex seating needs.

### Key features

- **Tilt-in-space and recline functions** – reduces shear and redistributes pressure (25° tilt, 30° recline)
- **Dual Stability System (DSS)** – enhances safety and balance during tilting
- **Customisable fit** – adjustable seat depth, backrest, armrest, headrest and legrest
- **Self-propelled and transit options** – accommodates different mobility needs
- **Durable frame with 135kg capacity** – supports long-term daily use
- **High manoeuvrability** – smooth indoor and outdoor movement, 77cm turning radius
- **Easy-clean PU stretch fabric** – supports hygiene and maintenance
- **Standard accessories** – includes headrest, seat belt, calf pads and lateral supports

## ADDITIONAL INFORMATION

<b>Brand</b>	Invacare
<b>Max SWL</b>	135 kg
<b>Seat depth</b>	43cm - 50cm
<b>Seat width</b>	37cm - 59cm depending on model
<b>Model</b>	14" Transit, 16" Self Propelled, 16" Transit, 18" Self Propelled, 18" Transit, 20" Self Propelled, 20" Transit
<b>Colour</b>	Pearl Grey