




**SKU:** S0055

## **PHILLIPS EVERFLO 5 L FLOOR OXYGEN CONCENTRATOR**

**\$2,699.99**

The Philips EverFlo 5 L Floor Oxygen Concentrator provides continuous oxygen therapy at flow rates up to 5 L/min, with oxygen purity at 93% ( $\pm$  3%). Its quiet, 14 kg design is easy to move and maintain, making it ideal for home based respiratory therapy.



Scan QR code to  
view product online

## PRODUCT DESCRIPTION

The Philips EverFlo 5 L Floor Oxygen Concentrator is designed for clients requiring continuous oxygen therapy at home.

Delivering flow rates from 0.5 L/min to 5 L/min, it maintains oxygen purity at 93% ( $\pm 3\%$ ), ensuring consistent medical-grade therapy. The integrated Oxygen Percentage Indicator (OPI) provides real-time monitoring of oxygen concentration, helping track performance easily.

Its compact 14 kg design, combined with castors and a recessed flow meter, supports safe repositioning and user-friendly operation. The device features intuitive controls, bright LEDs, and a durable metal cannula connector to enhance safety.

With its low-maintenance design, the EverFlo offers reliable, everyday therapy while reducing service requirements, making it a practical solution for long-term use.

### Who it's for

- For clients requiring continuous oxygen therapy at home, particularly those with chronic respiratory conditions such as COPD.
- For therapists seeking a reliable, low-maintenance oxygen concentrator for patients needing up to 5 L/min continuous flow.

### Key features

- **Continuous flow up to 5 L/min** – accommodates those with varying oxygen requirements
- **Consistent oxygen concentration** – maintains 93%  $\pm 3\%$  at 5 L/min
- **Compact and lightweight** – weighs only 14 kg, facilitating easy transport and placement
- **Integrated OPI technology** – ensures accurate oxygen delivery and enhances patient safety
- **User-friendly design** – features intuitive controls, bright LEDs, and a protected cannula fitting
- **Enhanced safety features** – includes a recessed flow meter and durable metal cannula connector

## ADDITIONAL INFORMATION

<b>Flow Range</b>	0-5 LPM
<b>Sound Level</b>	51.5 dBA
<b>Size (H x W x D)</b>	584 mm x 381 mm x 241 mm
<b>Weight</b>	14 kg
<b>Power Consumption</b>	590 Watt
<b>Brand</b>	Philips Respironics