



SKU: S0055

PHILLIPS EVERFLO 5 L FLOOR OXYGEN CONCENTRATOR

\$2,699.99

The Philips EverFlo 5 L Floor Oxygen Concentrator provides continuous oxygen therapy at flow rates up to 5 L/min, with oxygen purity at 93% ($\pm$ 3%). Its quiet, 14 kg design is easy to move and maintain, making it ideal for home based respiratory therapy.



Scan QR code to
view product online

PRODUCT DESCRIPTION

The Philips EverFlo 5 L Floor Oxygen Concentrator is designed for clients requiring continuous oxygen therapy at home.

Delivering flow rates from 0.5 L/min to 5 L/min, it maintains oxygen purity at 93% ($\pm 3\%$), ensuring consistent medical-grade therapy. The integrated Oxygen Percentage Indicator (OPI) provides real-time monitoring of oxygen concentration, helping track performance easily.

Its compact 14 kg design, combined with castors and a recessed flow meter, supports safe repositioning and user-friendly operation. The device features intuitive controls, bright LEDs, and a durable metal cannula connector to enhance safety.

With its low-maintenance design, the EverFlo offers reliable, everyday therapy while reducing service requirements, making it a practical solution for long-term use.

Who it's for

- For clients requiring continuous oxygen therapy at home, particularly those with chronic respiratory conditions such as COPD.
- For therapists seeking a reliable, low-maintenance oxygen concentrator for patients needing up to 5 L/min continuous flow.

Key features

- **Continuous flow up to 5 L/min** – accommodates those with varying oxygen requirements
- **Consistent oxygen concentration** – maintains 93% $\pm 3\%$ at 5 L/min
- **Compact and lightweight** – weighs only 14 kg, facilitating easy transport and placement
- **Integrated OPI technology** – ensures accurate oxygen delivery and enhances patient safety
- **User-friendly design** – features intuitive controls, bright LEDs, and a protected cannula fitting
- **Enhanced safety features** – includes a recessed flow meter and durable metal cannula connector

ADDITIONAL INFORMATION

Weight	14 kg
Length	38.1 cm
Width	24.1 cm
Height	58.4 cm
Oxygen Flow	Continuous Flow
Flow Rate	0.5 – 5 LPM
Noise Level	45 dBA*
Operation	Simple control functions
Brand	Philips Respironics