




RESPIRONICS EVERFLO OXYGEN CONCENTRATOR - 5 LITRE

~~\$1,594.99~~ **\$1,499.99**

Until now, oxygen concentrators have been very similar - heavy, bulky, noisy, or require frequent maintenance. The EverFlo from Respironics is a unique stationary concentrator that delivers the features home care providers want and patients deserve. Doctor's prescription may be required for oxygen concentrator.



SKU: S0055

Categories: [Floor Oxygen Concentrators](#),
[Oxygen Concentrator](#)

GALLERY IMAGES



PRODUCT DESCRIPTION

This Respironics Everflo unit is one of the quietest and lightest concentrators available in Australia. It will provide up to 5 Litres of Oxygen per minute. It is designed to be used in a domestic setting and operates on 240 V power.

Please note that a prescription is required for you to be able to purchase medical oxygen.

ADDITIONAL INFORMATION

| | |
|-----------------------------------|--------------------------|
| Flow Range | 0-5 LPM |
| Sound Level | 51.5 dBA |
| Size (H x W x D) | 584 mm x 381 mm x 241 mm |
| Weight | 14 kg |
| Power Consumption | 590 Watt |
| Brand | Philips Respironics |
| Manufacturer's part number | 1020010 |