



SKU:

## REVOLUTION SWIVEL LOWBACK DAY CHAIR

**\$1,249.99**

The Revolution Swivel Low Back Day Chair features a patented swivel slide mechanism and adjustable seat to support safer transfers, whether independent or assisted. Designed for comfort, functionality and hygiene, it's suited for both home and care environments.



Scan QR code to  
view product online

VARIATIONS

Image	SKU	Price	Description	Colour
	S14475	\$1,249.99		Sandstone
	S326264	\$1,249.99		Espresso
	S310793	\$1,249.99		Cloud

## PRODUCT DESCRIPTION

The Revolution Swivel Low Back Day Chair is a practical seating solution designed to improve safety and independence during transfers. It features a patented swivel slide mechanism that allows the seat to rotate 360° and slide forwards. This allows individuals to move closer to a table or desk safely, without standing.

Operated via a user-friendly paddle lever, the swivel mechanism supports smoother transitions in and out of the chair. The chair also offers two seat depth settings and four height adjustment positions to help achieve optimal postural alignment and fit.

Made from Terra Nova Neo commercial-grade fabric, the chair offers a soft leather-like feel while providing superior stain and wear resistance. The surface is easily cleaned, making it ideal for shared or clinical settings.

### Who it's for

- For clients seeking a static chair with swivel and sliding depth adjustability for easier transfers
- For therapists seeking a durable, functional day chair with a swivel-slide mechanism to support independent transfers

### Key features

- **Patented swivel slide mechanism** – improves access and transfer safety with 360° rotation and forward glide
- **Easy-to-use hand lever** – allows controlled swivel and slide movement
- **Height-adjustable legs** – four settings to customise seating height
- **Sliding seat depth** – seat slides forwards allowing closer access to tables or other objects
- **Terra Nova Neo fabric** – soft-touch finish with commercial-grade stain and wear resistance

## **ADDITIONAL INFORMATION**

<b>Overall Width</b>	505 mm
<b>Overall Depth</b>	560 mm
<b>Seat Height</b>	475 mm – 550 mm
<b>Seat Width</b>	445 mm
<b>Seat Depth</b>	480 mm
<b>Maximum User Weight</b>	120 kg