




**SKU:**

## **RIFTON ACTIVITY CHAIR**

The Rifton Activity Chair provides exceptional clinical positioning for children and adults.

The Activity Chair can be used for children with simple to complex needs, for everyday activities such as feeding, speech therapy and active learning. It is also suitable for children with sensory processing challenges.



Scan QR code to  
view product online

## **PRODUCT DESCRIPTION**

It is available with a Standard or Hi-Lo base. The Hi-Lo Base provides excellent positioning and comfort while maximising social inclusion and convenience.

The Rifton Activity chair features tool-free adjustments, and can be used for most home/classroom activities at desk or table height. The angles of the tilt in space and backrest recline can be adjusted while your child is in the chair.

The Activity Chair is available in three sizes, and each can accommodate a generous growth range. A large range of accessories is also available.

Available in three sizes

Max user weight 102 kg

### **Standard Base**

- Seat and back padding including adjustable back and seat belt
- Base with tilt-in-space
- Armrests and short adjustable legs
- Caster swivel lock
- Choice of 6 pad colors

### **Hi/Lo Base**

- Seat and back padding including adjustable back, seat belt and footboard
- Wheeled base with tilt-in-space and height adjustable
- Armrests
- Push handles

