


SKU:

## ROHO COMMUNE CUSHION

~~\$1,399.99~~ \$1,399.99



The Roho Commune Cushion is made up of interconnected air cells and has a unique design that promotes blood flow. The Commune Cushion has a section cut out in the centre of the cushion for toileting. The flexible cells will adapt to the contours of the user's body, reducing any shear forces while repositioning or transferring.



Scan QR code to  
view product online

VARIATIONS

Image	SKU	Price	Description	Size
	S501	\$1,124.99		15x16"
	S7724	\$1,399.99		16x16"
	S7725	\$1,124.99		16x18"
	S7726	\$1,124.99		18x18"

## **PRODUCT DESCRIPTION**

The Roho Commode Cushion is made up of durable neoprene air cells that provide a high quality of pressure care, making it excellent for those who suffer from skin breakages, or pressure sores. As well as being preventative these cushions can aid with the healing process of pressure sores.

- Cells are interconnected and their unique shape promote blood circulation to facilitate pressure prevention, healing and relief
- Ideal for clients who are at high risk of skin breakdown
- Air pressure is adjustable and can be customised to suit the user
- Durable neoprene construction is soft to suit skin sensitivity and is easily cleaned
- Flexible cells adapt to patient contours and reduce shear forces while transferring or re positioning

## **ADDITIONAL INFORMATION**

<b>Height</b>	60mm
<b>Max SWL</b>	Unlimited
<b>Brand</b>	Permobil