




SKU:

## **PRESSURE CUSHION - ROHO SINGLE VALVE LOW PROFILE**

**~~\$1,239.99~~ \$1,239.99**

The ROHO Low Profile Pressure Cushion features interconnected cells with a unique shape which promotes blood circulation to facilitate pressure prevention, healing and relief. Flexible cells adapt to patient contours and reduce shear forces while transferring or repositioning.



Scan QR code to  
view product online

**VARIATIONS**

Image	SKU	Price	Description	Dimensions
	S17001	\$1,099.99		22" x 22"
	S301139	\$1,239.99		24" x 18"
	S12350	\$1,099.99		24" x 20"
	S12348	\$1,079.99		13" x 13"
	S2413	\$1,079.99		15" x 15"
	S2300	\$1,079.99		15" x 16"
	S2094	\$1,079.99		16" x 16"
	S2871	\$1,079.99		16" x 18"
	S2870	\$1,079.99		16" x 20"
	S2090	\$1,079.99		18" x 16"
	S2684	\$1,079.99		18" x 18"
	S9891	\$1,079.99		18" x 20"
	S11740	\$1,079.99		20" x 16"
	S2092	\$1,079.99		20" x 18"
	S9858	\$1,079.99		20" x 20"
	S1672	\$1,239.99		22" x 18"

## **PRODUCT DESCRIPTION**

The cellular design and soft, malleable material allow the user to be immersed in the cushion for optimal pressure redistribution while reducing friction and shear. The easy adjustment enables the cushions to conform to the individual's shape, minimising tissue deformation and optimising blood circulation.

## ADDITIONAL INFORMATION

<b>Brand</b>	Permobil
<b>Height</b>	65 mm
<b>Max SWL</b>	Unlimited