



SKU:

ROHO MOSAIC CUSHION WITH HD COVER

\$399.99

The Roho Mosaic Cushion offers affordable skin protection and comfort for those at low risk of pressure injury. With interconnected air cells that adapt to each person's shape, it provides a lightweight, supportive seating option for everyday use.



Scan QR code to
view product online

VARIATIONS

Image	SKU	Price	Description	Sizes
	S12439	\$399.99		16" x 16"
	S12438	\$399.99		18" x 16"
	S12437	\$399.99		18" x 18"

PRODUCT DESCRIPTION

The Roho Mosaic Cushion is a cost-effective, air-filled seating solution designed to support skin health and comfort for individuals at low risk of pressure injury. Its segmented design with interconnected air cells conforms to the seated shape, helping to distribute pressure evenly and reduce the risk of skin breakdown.

Constructed from flame-retardant PVC, the Roho Mosaic Cushion includes a two-way stretch cover with a non-skid bottom to keep it securely in place during use.

The lightweight design ensures easy handling and transport, while the included hand pump allows quick and easy inflation and adjustment.

Who it's for

- For clients who need a lightweight, affordable cushion to help reduce pressure injury risk in daily use.
- For therapists prescribing a basic, adjustable cushion that promotes comfort and skin health for clients in community and home settings.

Key features

- **Interconnected air cells** – adapt to each person's shape to distribute pressure evenly
- **Adjustable air inflation** – allows customisation for individual comfort and support
- **Lightweight design** – easy to handle, transport and maintain
- **Flame-retardant PVC construction** – promotes durability and safety
- **Heavy-duty cover** – includes non-skid bottom for secure placement
- **Hand pump and repair kit included** – supports easy setup and ongoing maintenance

ADDITIONAL INFORMATION

Sizes

16" x 16", 18" x 16", 18" x 18"