



SKU:

ROYALE PRESSURE RELIEVING LIFT CHAIR (TWIN MOTOR)




\$3,199.99

Exclusive to ILS, the Royale Air Pressure Relieving Lift Chair is designed to accommodate extended seating periods and is intended for use by individuals at risk of pressure injury. Featuring a unique air bladder system within the headrest, lumbar and seat cushions, the chair is upholstered with a medical-grade fabric designed to help reduce the risk of pressure injuries. The Royale Air Pressure Relieving Lift Chair has been clinically tested to reduce pressure in common seating positions. The chair is easily moved and positioned with front and rear castors.



Scan QR code to
view product online

VARIATIONS

Image	SKU	Price	Description	Sizes
	S307655	\$3,199.99		Small
	S307656	\$3,199.99		Medium
	S307657	\$3,199.99		Large

PRODUCT DESCRIPTION

The Royale Air Pressure Relieving Lift Chair is purpose-built for individuals who are seated for extended periods and at risk of developing pressure injuries. Its unique air bladder system – integrated into the headrest, lumbar and seat cushions – can be adjusted via discreet hand pumps to meet individual pressure relief needs.

Upholstered in soft, breathable four-way stretch medical-grade fabric, the chair helps protect skin integrity across all key contact areas.

The twin motor mechanism provides independent control of the backrest and legrest, to accommodate both individual comfort and clinical postural needs.

Swivel castors with rear brakes allow easy chair movement and secure placement once positioned. Additional features include a backlit handset with USB charging and accessible side pockets.

Who it's for

- For clients who require pressure management support during long periods of seated care.
- For therapists needing a clinically tested lift chair with customisable positioning and mobility features.

Key features

- **Adjustable air bladder system** – supports headrest, lumbar and seat pressure relief
- **Hand pumps** – allow pressure adjustment and store neatly when not in use
- **Medical-grade stretch fabric** – breathable and soft to reduce skin shear
- **Twin motor mechanism** – independent control of backrest and legrest
- **Swivel castors with rear brakes** – allow for easy positioning
- **Backlit hand control with USB** – simple to use and convenient for charging devices
- **Dual side storage pockets** – keep essentials close at hand
- **Battery back-up system** – ensures continuous operation during power loss

ADDITIONAL INFORMATION

Brand	Royale Medical
Max SWL	160 kg
Backrest Height	Small: 680 mm / Medium: 720 mm / Large: 740 mm
Seat Height	Small: 480 mm / Medium: 520 mm / Large: 520 mm
Seat Depth	Small: 500 mm / Medium: 500 mm / Large: 530 mm
Seat Width	Small: 500 mm / Medium: 500 mm / Large: 595 mm
Overall Width	Small: 750 mm / Medium: 750 mm / Large: 820 mm
Overall Height	Small: 1060 mm / Medium: 1060 mm / Large: 1060 mm
Sizes	Small, Medium, Large
Warranty	2 years