



SKU: S19567

ROYALE CHADWICK OXFORD PLUSH LEATHER LIFT CHAIR (QUAD MOTOR) - MARINE BLUE

\$4,399.99

Exclusive to ILS, the Royale Chadwick combines advanced quad motor technology with luxurious 'Oxford Plush' leather upholstery, offering the ultimate seating solution. The Royale Chadwick has fully independent powered headrest, lumbar support, backrest and legrest functions to provide users with the greatest level of control and flexibility in positioning needs.



Scan QR code to
view product online

PRODUCT DESCRIPTION

The Royale Chadwick Oxford Plush Leather Lift Chair (Quad Motor) provides both comfort and style. Its quad motor system allows full independent control of the headrest, lumbar, backrest and legrest, while the Oxford Plush leather upholstery finish offers a softer feel for those preferring a more luxurious finish.

Crafted from high-quality hides, the Oxford Plush Leather is breathable, durable and soft to the touch. The chair also includes ergonomic features such as a two-section backrest, multi-layered foam and a backlit control with USB port.

Who it's for

- For clients who want a softer-feel leather chair with premium finish and full-body adjustability.
- For therapists looking for a high-spec quad motor chair that supports comfortable posture and simple adjustment for their client.

Key features

- **Quad motor adjustment** – independently controls headrest, lumbar, backrest and legrest
- **Oxford Plush Italian leather** – soft, breathable and durable
- **Two-section backrest** – supports comfort and alignment
- **Premium layered foam** – enhances long-term seating comfort
- **Handset with backlight and USB** – includes 'Home' return function
- **Side storage pockets** – keep essentials close by
- **Battery backup** – maintains functionality during power loss

ADDITIONAL INFORMATION

Brand	Royale Medical
Max SWL	150 kg (330.6 lbs)
Overall height	110 cm (43.3 inches)
Overall width	85 cm (33.4 inches)
Overall depth	80 cm (31.4 inches)
Armrest height (from floor)	64 cm (25 inches)
Seat height (from floor)	50 cm (19.6 inches)
Seat width (measured on sitting area)	56 cm (22 inches)
Seat depth (measured on sitting area)	50 cm (19.6 inches)
Max seat height (lifted)	57 cm (22.4 inches)
Backrest length (fully reclined)	60 cm (23.6 inches)
Max legrest height (lifted)	52 cm (20.4 inches)
Foam density	35/70
Upholstery	Genuine leather
Colour options	Colonial hemp beige, dark cigar brown, jet carbon black, marine blue
Warranty	Electrics including transformer, hand control, motor, and all cables – 2 years Steel lift/recline mechanism – 2 years