



ROYALE MEDICAL CHADWICK SOFT TOUCH FABRIC LIFT CHAIR

\$2,395.99

Exclusive to the ILS group, the state-of-the-art Royale Medical Chadwick recliner lift chair puts you in complete control of your comfort. Its unique quad motor technology ensures fully independent headrest, lumbar, backrest and footrest adjustments.

Boasting a contemporary and elegant design, the Chadwick recliner lift chair easily matches your living room furniture no matter your decor style. It is upholstered in premium soft touch fabric and available in three colours - indigo blue, oyster latte and coffee bean brown. Experience ultimate luxury and ergonomic satisfaction every day with the Chadwick lift chair.



SKU: N/A

Categories: [Lift Chairs](#), [Quad Motor Lift Chairs](#)

VARIATIONS

| Image | SKU | Price | Description | Colour |
|---|--------|------------|-------------|-------------------|
|  | S19401 | \$2,395.99 | | Indigo blue |
|  | S19400 | \$2,395.99 | | Oyster latte |
|  | S19405 | \$2,395.99 | | Coffee bean brown |

GALLERY IMAGES



PRODUCT DESCRIPTION

The Royale Medical Chadwick recliner lift chair offers ultimate support, style, durability and comfort.

Key features include:

- Quad motors - Adjust the headrest, lumbar, backrest and footrest independently to suit your positioning needs.
- Multilayered foam - The highest quality multilayered foam is used for ultimate comfort and product longevity.
- Soft touch fabric upholstery in three colour options of indigo blue, oyster latte and coffee bean

brown.

- Home button – Press the unique home button feature to bring the chair back to its neutral position, eliminating the hassle of multiple operations.
 - USB – Charge your electrical devices through the USB ported hand control.
 - Soft vertical lift – Move safely and easily out of the Chadwick recliner lift chair.
- Transport castors – Move the Chadwick recliner lift chair with ease using the transport castors provided.
- Warranty period of 2 years for electrical and steel lift/recliner mechanisms and 5 years for the timber frame.

ADDITIONAL INFORMATION

| | |
|--|----------------------|
| Brand | Royale Medical |
| Max SWL | 150 kg (330.6 lbs) |
| Overall height | 110 cm (43.3 inches) |
| Overall width | 85 cm (33.4 inches) |
| Overall depth | 80 cm (31.4 inches) |
| Armrest height (from floor) | 64 cm (25 inches) |
| Seat height (from floor) | 50 cm (19.6 inches) |
| Seat width (measured on sitting area) | 56 cm (22 inches) |
| Seat depth (measured on sitting area) | 50 cm (19.6 inches) |
| Max seat height (lifted) | 57 cm (22.4 inches) |
| Backrest length (fully reclined) | 60 cm (23.6 inches) |

| | |
|------------------------------------|---|
| Max legrest height (lifted) | 52 cm (20.4 inches) |
| Foam density | 35/70 |
| Upholstery | Soft touch fabric |
| Colour options | Indigo blue, oyster latte, coffee bean brown |
| Warranty | Electrics including transformer, hand control, motor, and all cables - 2 years Steel lift/recline mechanism - 2 years Timber frame (structural and workmanship) - 5 years |