




**SKU:**

## **ROYALE MAYFAIR MINI LUXURY LEATHER (TWIN MOTOR)**

**\$0.00**

Exclusive to ILS, the Royale Mayfair Mini Luxury leather lift chair combines premium materials with quality twin motor functionality, offering individuals a comfortable, supportive and attractive lift chair solution. Compact in size for smaller users, the Mayfair Mini Luxury is designed with a twin motor system to provide users independent control of the backrest and legrest functions to ensure flexible and comfortable positioning.



Scan QR code to  
view product online

## VARIATIONS

**Image**

**SKU**

**Price**

**Description**

**Colors**



S300955 \$2,699.99

Espresso Dark Brown



S300956 \$2,699.99

Taupe Beige

## PRODUCT DESCRIPTION

The Royale Mayfair Mini Luxury Leather Lift Chair (Twin Motor) is designed for smaller individuals seeking additional support and assistive movement. Its twin motor system allows the backrest and legrest to be positioned independently for flexible comfort.

Upholstered in breathable, high-quality leather from premium hides, the chair also features a supportive three-section backrest and multi-layered foam.

Its compact dimensions make it well suited to petite individuals and smaller living spaces.

### Who it's for

- For clients under 167 cm who want a luxurious, compact lift chair with flexible positioning.
- For therapists supporting petite clients who need a twin motor lift chair with leather upholstery and strong ergonomic support.

### Key features

- **Twin motor adjustability** – backrest and legrest can be moved independently
- **Compact sizing** – ideal for individuals 167 cm and under
- **Luxury leather finish** – soft, durable and breathable
- **Three-section waterfall backrest** – supports spinal alignment and comfort
- **Quality foam seat** – provides pressure relief and structure
- **Battery back-up** – maintains functionality during outages
- **Intuitive handset and side pockets** – for everyday convenience
- **Side pockets** – practical and accessible storage

## ADDITIONAL INFORMATION

<b>Brand</b>	Royal Medical
<b>Max SWL</b>	150 kg (330.6 lbs)
<b>Motor</b>	Twin
<b>Overall height</b>	98 cm (38.5 inches)
<b>Overall width</b>	66 cm (25.9 inches)
<b>Overall depth</b>	75 cm (29.5 inches)
<b>Armrest height (from floor)</b>	56 cm (22 inches)
<b>Seat height (from floor)</b>	45 cm (17.7 inches)
<b>Seat width (measured on sitting area)</b>	46 cm (18.1 inches)
<b>Seat depth (measured on sitting area)</b>	43 cm (16.9 inches)
<b>Seat max height (lifted)</b>	57 cm (22.4 inches)
<b>Back length (fully reclined)</b>	19 cm (7.4 inches)
<b>Max legrest height (lifted)</b>	47 cm (18.5 inches)
<b>Foam density</b>	28/70
<b>Upholstery</b>	Premium bovine leather
<b>Colour options</b>	Espresso dark brown and taupe beige
<b>Warranty</b>	2 years for electrical and steel lift/recliner mechanisms and frame integrity.