



**SKU:**

## **ROYALE MAYFAIR MINI LUXURY LEATHER (TWIN MOTOR)**


**\$0.00**

Designed to be small, the Royale Medical's Mayfair Mini Luxury lift chair is compact in size for petite people. It is packed with quality features much loved by customers of the Mayfair luxury lift chair such as a quiet and smooth three position twin motor lift system and very comfortable 3-section *waterfall* backrest design for optimum comfort and back support. Available in two colours, it is made with premium bovine leather to suit any setting or home decor. If you're a petite person, the Mayfair Mini Luxury recliner lift chair provides customised comfort and support specifically designed for your size!



Scan QR code to  
view product online

### VARIATIONS

| Image  | SKU     | Price      | Description | Colors              |
|--|---------|------------|-------------|---------------------|
|   | S300955 | \$2,699.99 |             | Espresso Dark Brown |
|  | S300956 | \$2,699.99 |             | Taupe Beige         |

## PRODUCT DESCRIPTION

Royale Medical's Mayfair Mini Luxury recliner lift chair offers support, comfort and style for every day use at an affordable price. The compact size has been specifically designed for petite people and is also fantastic for small rooms.

Key features include:

- Very comfortable 3-section *waterfall* backrest design for optimum comfort and upper, middle and lower back support.
- Easy-to-use functions so you can easily recline the backrest and legrest.
- Quiet and smooth twin motor lift system.
- Quality actuators and American-made mechanisms allows motors to run quietly for a smooth and discreet operation.
- Quality foams and no sag springs providing you with supreme postural support.
- Premium bovine leather upholstery in two colours - Choose either espresso dark brown and taupe beige.
- Backup power supply for a peace of mind.
- Side storage pocket.
- Manufacturer's warranty period of 2 years for electrical and steel lift/recliner mechanisms and frame integrity.

### Single vs twin vs quad motor lift chairs - what's the difference?

All ILS lift chairs have a lifting action, to support the user to stand from a seated position.

- Single motor chairs move the back rest and leg rest at the same time.
- Twin motor chairs allow the backrest and leg rest to move independently of each other, typically allowing for a greater degree of leg elevation and greater degree of recline compared to a standard single motor chair.
- Quad motor chairs allow independent movement of headrest, lumbar, backrest and footrest allowing for the greatest level of re-positioning of all chairs.

## ADDITIONAL INFORMATION

|  |  |
|--|--|
| <b>Brand</b>                                 | Royal Medical  |
| <b>Max SWL</b>                               | 150 kg (330.6 lbs)   |
| <b>Motor</b>                                 | Twin   |
| <b>Overall height</b>                        | 98 cm (38.5 inches)  |
| <b>Overall width</b>                         | 66 cm (25.9 inches)  |
| <b>Overall depth</b>                         | 75 cm (29.5 inches)  |
| <b>Armrest height (from floor)</b>           | 56 cm (22 inches)  |
| <b>Seat height (from floor)</b>              | 45 cm (17.7 inches)  |
| <b>Seat width (measured on sitting area)</b> | 46 cm (18.1 inches)  |
| <b>Seat depth (measured on sitting area)</b> | 43 cm (16.9 inches)  |
| <b>Seat max height (lifted)</b>              | 57 cm (22.4 inches)  |
| <b>Back length (fully reclined)</b>          | 19 cm (7.4 inches)   |
| <b>Max legrest height (lifted)</b>           | 47 cm (18.5 inches)  |
| <b>Foam density</b>                          | 28/70  |
| <b>Upholstery</b>                            | Premium bovine leather   |
| <b>Colour options</b>                        | Espresso dark brown and taupe beige  |
| <b>Warranty</b>                              | 2 years for electrical and steel lift/recliner mechanisms and frame integrity. |