



SKU:

ROYALE MEDICAL CHADWICK MINI LIFT CHAIR

\$0.00

Exclusive to Independent Living Specialists, the Chadwick Collection features unique quad motor technology, with fully integrated, independent powered headrest, lumbar, backrest and footrest adjustments, putting you in complete control of your comfort and positioning needs.

Boasting a contemporary and elegant design, the Chadwick Collection suits all décor styles and was designed with comfort and ergonomics in mind.

Recommended for individuals 167cm/5'6" and below.



Scan QR code to
view product online

VARIATIONS

Image	SKU	Price	Description	Type of material	Color
	S306218	\$3,999.99		Oxford Plush Leather	Oxford Plush Leather Marine Blue
	S306213	\$3,999.99		Oxford Plush Leather	Oxford Plush Leather Jet Carbon Black
	S306215	\$3,399.99		Genuine Leather	Leather Espresso Dark Brown
	S306216	\$3,399.99		Genuine Leather	Leather Taupe Beige
	S306217	\$2,899.99		Soft Touch Fabric	Soft Coffeebean Brown
	S306214	\$2,899.99		Soft Touch Fabric	Soft Indigo Blue

PRODUCT DESCRIPTION

Key features

- Soft vertical lift action to enable safe, easy movement and transfer out of the chair.
- Quad motor allowing for fully independent headrest, lumbar, backrest and footrest adjustments.
- 2-year warranty on all electrical components and 5-year structural warranty (See specifications for more information).
- Choice of 3 premium materials, each with colour options.
- Highest quality multi-layered foam for ultimate comfort and durability.
- Backup power supply for peace of mind.
- Hand control with unique home button bringing the chair back to its neutral position in one operation.
- USB port on the hand control to charge your electrical devices.

Material

Soft Touch Fabric – This scratch resistant material is breathable, warming and soft for comfort.

Genuine Leather – This natural material has hypo-allergenic properties. It is durable, breathable, easy to maintain and attractive.

Oxford Plush Leather – With all the benefits of Genuine Leather, this Italian leather is made with high quality hides producing a softer finish for increased comfort.

Single vs twin vs quad motor lift chairs - what's the difference?

All ILS lift chairs have a lifting action, to support the user to stand from a seated position.

- Single motor chairs move the back rest and leg rest at the same time.
- Twin motor chairs allow the backrest and leg rest to move independently of each other, typically allowing for a greater degree of leg elevation and greater degree of recline compared to a standard single motor chair.
- Quad motor chairs allow independent movement of headrest, lumbar, backrest and footrest allowing for the greatest level of re-positioning of all chairs.

ADDITIONAL INFORMATION

Brand	Royale Medical
Max SWL	150 kg (330.6 lbs)
Overall height	102 cm (40.2 inches)
Overall width	75 cm (29.5 inches)
Overall depth	72 cm (28.3 inches)
Armrest height (from floor)	60.5 cm (23.8 inches)
Seat height (from floor)	44 cm (17.3 inches)
Seat width (measured on sitting area)	45 cm (17.7 inches)
Seat depth (measured on sitting area)	51 cm (20.1 inches)
Max seat height (lifted)	57 cm (22.4 inches)
Backrest length (fully reclined)	60 cm (23.6 inches)
Max legrest height (lifted)	52 cm (20.4 inches)
Foam density	35/70
Upholstery	Genuine leather , Soft Touch Fabric , Oxford Plush Leather
Colour options	Taupe Beige, Espresso Dark Brown, Indigo Blue, Coffeebean Brown, Marine Blue, Jet Cabon Black
Warranty	Electrics including transformer, hand control, motor, and all cables and steel lift/recline mechanism – 2 years Timber frame – 5 years