


SKU:

SUNRISE MEDICAL BREEZY RUBIX WHEELCHAIR

\$1,999.99




The Sunrise Medical Breezy Rubix Wheelchair is a folding, heavy-duty wheelchair that supports individuals up to 170 kg. With seat widths up to 24" and multiple adjustable features, it offers tailored comfort, postural support and reliable daily mobility.



Scan QR code to
view product online



VARIATIONS

Image	SKU	Price	Description	Sizes
	S10000	\$1,999.99		20"
	S10001	\$1,999.99		22"
	S10002	\$1,999.99		24"

PRODUCT DESCRIPTION

The Breezy Rubix Wheelchair is designed for individuals needing wider seating and a higher weight capacity. It features a reinforced aluminium frame, adjustable seat depth, configurable seat height, and a tension-adjustable backrest for tailored support.

Quick-release rear wheels and swing-away, angle-adjustable legrests support easier transfers and compact storage.

Additional features include flip-back armrests, anti-tip wheels and rear wheel park brakes for everyday usability and safety. It also has crash test certification.

Who it's for

- For clients seeking a reliable, adjustable heavy-duty wheelchair with wider seat width and self-propulsion.
- For therapists prescribing a bariatric wheelchair with certified safety, daily durability and adjustability.

Key features

- **Adjustable seat depth (46 cm - 51 cm)** - supports postural alignment and comfort
- **Heavy-duty aluminium frame** - accommodates up to 170 kg
- **Available seat widths** - 20", 22" and 24" accommodates larger body shapes
- **Flip-back armrests** - allow safer lateral transfers and access
- **Tension-adjustable backrest** - offers personalised support
- **Swing-away, angle-adjustable legrests** - enhance positioning and transfers
- **Quick-release rear wheels** - for ease of transport and handling
- **Rear park brake** - accessible and easy to operate
- **Crash tested** - certified for transport safety

ADDITIONAL INFORMATION

Brand	Sunrise Medical
Seat angle	0° to 6°
Seat tilt	No
Product weight	From 16.8 Kg
Sizes	20", 22", 24"